Myth Busting

AAC in Early Intervention

AbleNet

QuickTalker Freestyle[®]
LET THE CONVERSATION BEGIN

1

Myth: "My child will become dependent on their AAC device and lose their voice"

FACT: Research shows AAC supports and often enhances spoken language by giving children a reliable way to communicate.

2

Myth: "AAC devices are too high-tech for young children"

FACT: Devices can be customized to start simple and grow with your child's communication skills.

3

Myth: "You have to be a certain age to get started with AAC"

FACT: Early use leads to better outcomes. AAC can support in developing language, literacy, and cognitive skills.

4

Myth: "Children this young shouldn't be engaging in screen time"

FACT: AAC devices are not the same as passive screen time. A speech-generating device should be treated as the child's voice, not for entertainment.

5

Myth: "There are prerequisites children need prior to using AAC"

FACT: Children do not need specific skills to begin AAC. Introducing AAC early supports language and cognitive development through early intervention.

Criscuola, M. (n.d.). Critical questions about screens, AAC, and kids with autism. ASHA Wire. Davidoff, B. (2017, January 1). AAC with energy—earlier - the ASHA leader. ASHA Wire. MaryAnn Romski, Rose A. Sevcik, Andrea Barton-Hulsey & Ani S. Whitmore (2015) Early intervention and AAC: What a difference 30 years makes. Augmentative and alternative communication (Baltimore, Md. : 1985). Romski, MaryAnn PhD, CCC-SLP; Sevcik, Rose A. PhD. Augmentative Communication and Early Intervention: Myths and Realities. Infants & Young Children 18(3):p 174-185, July 2005.

Aided Language Stimulation:

Supporting Early Intervention and Building AAC Skills Through Modeling

One of the best ways to support language in early intervention is to create a language rich environment, and it's no different for children using AAC! Research shows that the earlier children are exposed to their communication system, the better their long-term language outcomes (*Romski et al., 2010*). Think about how children learn to talk. They hear thousands of words before they say their first word. Children using AAC need the same rich exposure and consistent models to their AAC system from an early age.

What is Aided Language Stimulation?

Aided Language Stimulation (ALS), often referred to as modeling, is an evidence-based strategy where communication partners model language using an AAC system while speaking naturally (*Goossens, 1989*).

The Aided Language Stimulation technique:

1. Speak the complete sentence naturally - don't change how you normally talk
2. Point to 1-2 key words on the AAC device while you speak
3. No pressure or expectations - don't require the child to respond or imitate

“ Why Aided Language Stimulation Works:

Research demonstrates children who receive early, consistent aided language stimulation develop **stronger receptive language skills, better understanding of communication functions, and more robust expressive language abilities** over time (*Drager et al., 2006; Binger & Light, 2007*).

Snack

“I want more snack”
(Point to WANT + MORE)

Caregiving

“Let's put on shoes”
(Point to PUT + SHOES)

Transitions

“Time to go home”
(Point to GO + HOME)

Play

“My turn to build!”
(Point to MY + TURN)

Reading

“The bear is big”
(Point to BIG + BEAR)


Speech Sound Development

Speech sounds should be developed by the following ages. Speech development may be delayed if milestones are not met. There is a range for speech sound development and these sounds may be developed earlier than the provided age.

AbleNet

QuickTalker Freestyle
LET THE CONVERSATION BEGIN

Imitation Hierarchy + AAC Tips

Imitation is a core mechanism for language acquisition. This hierarchy moves from foundational motor imitation toward spoken and symbolic imitation. Progress is not always linear, and individuals may demonstrate skills at multiple levels simultaneously or show scatter across the hierarchy. AAC is a tool that can be used at every level!

1

Actions with Objects

Banging, pushing, stacking, knocking over

- ✓ Turn sound on — auditory feedback is powerful at this level.
- ✓ Bang the device like a drum to get attention, then model “boom!”
- ✓ Push the device across the table like a toy car, press “go!”
- ✓ Knock blocks over together, press “crash!” or “uh-oh!” in the moment
- ✓ Match sound words to the action in real time. Pairing the two is key!

2

Gross Motor Imitation

Clapping, stomping, arms up, patting knees

- ✓ Clap your hands and press “yay!” on the device at the same time.
- ✓ Stomp your feet and press “boom!” Pair the whole body movement with the word.
- ✓ Hold the device up and press “up!” while raising your arms together
- ✓ Keep it big, silly, and physical. Model the movement and the word simultaneously.

3

Facial & Oral Motor

Open mouth, tongue out, puff cheeks, blow

- ✓ Make a silly face and press “silly!” or “again!” right after.
- ✓ Blow out a candle or bubbles together, press “blow!” at the same moment.
- ✓ Use mirror play. Make the face side by side, then tap the device: “my turn,” “your turn”.

4

Gestures

Waving, pointing, nodding, shrugging

- ✓ Wave hello and press “hi!” simultaneously. Do both every time, together.
- ✓ When they point at something, immediately name it on the device
- ✓ Shrug and press “I don’t know!” Model communicative gestures as full messages.

5

Vocalizations & Sound Words

Animal sounds, “ooh/ahh,” sound effects

- ✓ Program sound words onto the device: “moo,” “crash,” “uh-oh,” “whoa,” “eww,” “yuck”.
- ✓ Make the sound out loud and press it on the device at the same time.
- ✓ Sound words are often the first words imitated. Prioritize and celebrate these sounds!
- ✓ Match the energy: big reaction, press the word, pause and wait.

6

Verbal Routines

“Ready, set... go!”, “1, 2... 3!”, song fill-ins

- ✓ Say the routine out loud and navigate to the final word on the device together.
- ✓ Pause before the last word, wait, then model on the device if needed.
- ✓ Keep energy and motivation high. Use these often and across preferred settings.

7

Single Words

Core vocabulary: more, help, stop, go, want

- ✓ Say the word out loud while pressing it, giving both models.
- ✓ Model in the moment it’s meaningful, not in isolation.
- ✓ Celebrate any approximation, whether spoken, on the device, or both!

8

Phrases

“want more,” “my turn,” “help me,” “all done”

- ✓ Say the phrase while navigating to each word on the device. Slow it down so they can follow.
- ✓ “Want,” press want, “more,” press more. Natural pace, together.
- ✓ Navigating two words is a motor planning win even before spontaneous use emerges.

Building Pre-Speaking Skills with AAC

Early communication starts long before first words. Skills like **turn-taking, imitation, and joint attention** can be the building blocks of language. High-tech AAC (Augmentative and Alternative Communication) can support these skills from the start. It's not a replacement for speech, but an additional way for individuals to connect and be understood.

1

Responding to the Environment

Noticing and Reacting to the World

- ✓ Turn sound on and **narrate the environment** in real time: phone rings, press "phone!"
- ✓ Comment on what's happening around you, not just what the individual is doing.
- ✓ Keep it **low pressure**. No response needed. You're building the association.

2

Responding to People

Showing awareness of and interest in communication partners

- ✓ Press "**hi!**" every time you greet them. Every time.
- ✓ Hold the device near your face so attention to you and the device happen together.
- ✓ **Wait after modeling**. Any response counts.

3

Turn-Taking

Learning that interaction goes back and forth

- ✓ Press a word, then **pause and wait**. The pause is the turn.
- ✓ Use fill-in-the-blank routines: "**Ready, set...**" press "**go!**"
- ✓ During a preferred song or game, press a word and leave space for any response.

4

Sustained Attention

Maintaining focus long enough for meaningful interaction

- ✓ Use AAC during whatever already holds their attention. **Follow their lead**.
- ✓ Short and successful beats long and effortful.
- ✓ End on a good moment. Leave them wanting more.

5

Joint Attention

Shared focus between a person and a referent

- ✓ When they look at something, immediately name it on the device.
- ✓ Press "**look!**" or "**wow!**" to draw shared attention to something interesting.
- ✓ Joint attention doesn't require eye contact. Follow their gaze, not just their eyes.

6

Purposeful Communication

Communicating with clear intent directed at someone

- ✓ Respond to every communication attempt as meaningful, regardless of how it looks.
- ✓ When they reach or point, model **“want ___!”** on the device right away.
- ✓ Pair your response with the device every time so they learn it works.

7

Following Directions

Demonstrating receptive understanding of language

- ✓ Press **“your turn”** before offering an opportunity.
- ✓ Press **“let’s go”** before every transition. Consistency builds comprehension.
- ✓ They likely understand far more than they can express. Keep that in mind.

8

Imitating Actions

Using imitation as a mechanism for learning language

- ✓ Model an action and press the word at the same time: knock blocks, press **“boom!”**
- ✓ Wave and press **“bye”** together. Always pair the motor with the language.
- ✓ Imitation can be immediate or delayed. Both count.
- ✓ See Imitation Hierarchy + AAC Tips for more information.

9

Initiating Communication

Leading the interaction rather than only responding

- ✓ Create genuine opportunities: place a preferred item just out of reach and wait.
- ✓ Wait with expectant attention before modeling or prompting.
- ✓ When they reach or signal, model **“I want ___”** or **“help please”** on the device.

Affirming Considerations for Autistic and Neurodiverse Individuals

- ✓ Sensory differences shape how someone engages. Work with them, not around them.
- ✓ Joint attention without eye contact is real and valid.
- ✓ Echolalia and scripts are meaningful communication. Respond to the intent.
- ✓ Motor planning differences may affect device access and imitation. Adapt the strategy, not the expectation.
- ✓ Assume competence. Model rich language regardless of current output.

Key Takeaways

- ✓ AAC supports communication at any age and any stage.
- ✓ All modalities are valid: gesture, vocalization, expression, device use.
- ✓ Model without pressure or expectation of immediate response.
- ✓ Follow the individual’s lead.
- ✓ Small, consistent progress is still progress.

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