




*Scope and Application:* For finishing baths. The estimated detection limit for program number 2350 is 0.05 g/L.



**1.** Press the soft key under **HACH PROGRAM**.

Select the stored program for autocatalytic nickel by pressing **2350** with the numeric keys.

Press: **ENTER**

**Note:** The Flow Cell and Sipper Modules cannot be used for this procedure.

**Note:** This method gives accurate results on most bath formulations. If the bath formulation in use responds differently, perform a new manual calibration. Prepare and store the calibration as directed under Calibration Standard Preparation.

**Note:** If samples cannot be analyzed immediately, see Sample Collection, Storage and Preservation following these steps.




**2.** The display will show: **HACH PROGRAM: 2350 Nickel, Autocatal.**

The wavelength ( $\lambda$ ), **720 nm**, is automatically selected.



**3.** Fill a sample cell with 10 mL of deionized water (the blank).



**4.** Fill a second cell with 10 mL of bath sample.


**Note:** Filter highly turbid samples.

**Note:** For proof of accuracy, use a 1000-mg/L (1.00-g/L) Nickel Standard Solution (see the OPTIONAL REAGENTS AND STANDARDS section) in place of the sample.




**5.** Add the contents of one Potassium 1 Reagent Powder Pillow to the sample (the prepared sample). Stopper and shake to dissolve.

**Note:** If a visible turbidity forms upon addition of Potassium 1 Reagent, dilute sample 1:1 with deionized water and repeat Step 5. Multiply results obtained in Step 8 by 2.




**6.** Place the blank into the cell holder. Close the light shield.



**7.** Press the soft key under **ZERO**.  
The display will show:

**0.00 g/L Ni<sup>2+</sup>**

**Note:** For alternate concentration units, press the **OPTIONS** soft key. Then press the soft key under **UNITS** to scroll through the available options. Press **ENTER** to return to the read screen.



**8.** Place the prepared sample into the cell holder. Close the light shield. The result in g/L Nickel (or chosen units) will be displayed.

## Interferences

**Table 1 Interfering Substances That Cause a Negative Interference**

Interfering Substance	Interference Levels and Treatments
Copper	All levels. Gives a similar blue color.

## Sample Collection, Storage and Preservation

Collect samples in clean plastic or glass bottles. Store at 4 °C (39 °F) or lower.

## Accuracy Check

### Standard Solutions Method

Use a 1000 mg/L Nickel Standard Solution in Step 4 in place of the sample. Perform the autocatalytic procedure as described above.

To adjust the calibration curve using the reading obtained with the 1000-mg/L Nickel Standard Solution, press the soft keys under **OPTIONS, MORE** then **STD:OFF**. Press **ENTER** to accept the displayed concentration, the value of which depends on the selected units. If an alternate concentration is used, enter the actual concentration and press **ENTER** to return to the read screen. See Section 1.5.5 *Adjusting the Standard Curve* for more information.

## Method Performance

### Precision

Standard: 1.00 g/L Ni<sup>2+</sup>

Program	95% Confidence Limits
2350	0.97–1.03 g/L Ni <sup>2+</sup>

For more information on determining precision data and method detection limits, refer to Section 1.5.

### Estimated Detection Limit

Program	EDL
2350	0.05 g/L Ni <sup>2+</sup>

For more information on derivation and use of Hach's estimated detection limit, see Section 1.5.2. To determine a method detection limit (MDL) as defined by the 40 CFR part 136, appendix B, see Section 1.5.1.

### Sensitivity

Program Number: 2350

Portion of Curve	ΔAbs	ΔConcentration
Entire Range	0.010	0.12 g/L Ni <sup>2+</sup>

See Section 1.5.3 *Sensitivity Explained* for more information.

## Calibration Standard Preparation

This method gives accurate results on most bath formulations. If the bath formulation in use responds differently, perform a new manual calibration.

- a. Prepare a sample of the finishing bath of known correct strength (100%). Or, titrate a sample of the bath to determine its exact strength.
- b. Dilute this bath sample 1:1 with deionized water to make a half-strength standard.
- c. Into three separate cells, add 10-mL of the following: deionized water (reagent blank), half-strength standard, and full-strength standard (100%).
- d. Using the colorimetric method for autocatalytic nickel and the calibration procedure described in the *User-Entered Programs* section of the *DR/4000 Spectrophotometer Instrument Manual*, generate a calibration curve from the standards prepared. Use the following user-method selections:

**Method name:** Nickel, Autocatal.

**Format:** XX.XX

**Units:** g/L

**Chemical form:** Ni

**Wavelength:** 720 nm

- e. Scroll to the Calibration table selection and press the soft key under **CREATE TABLE**. If the program has been stored previously, the soft key will display **EDIT TABLE** instead.

# NICKEL, AUTOCATALYTIC, continued

- f. Add the contents of one Potassium 1 Reagent Powder Pillow to each of the three cells as described in Step 5.
- g. Enter the nickel concentration of the reagent blank. Zero the instrument using deionized water contained in a sample cell. Measure the absorbance of the reagent blank and accept the reading by pressing the soft key under **READ**. Repeat this procedure for both the half and full-strength standards.

**Note:** Some variations of the calibration procedure are possible. See the DR/4000 Instrument Manual for details.

- h. When repeating the calibration for new bath formulations, select the same user-assigned program number in the procedure above. Prepare a new calibration curve by editing only the calibration table each time. Use the same stored program number and method settings.

This procedure provides good readability on baths up to 8 g/L. If the bath in use is more concentrated, dilute the sample and apply an appropriate correction factor.

## Summary of Method

A strong complexing agent chelates the nickel ions present in an “electroless” nickel bath to form a blue colored chelate. The blue color is then measured directly to give the g/L nickel present in the bath.

## Safety

Good safety habits and laboratory techniques should be used throughout the procedure. Consult the *Material Safety Data Sheet* for information specific to the reagents used. For additional information, refer to Section 3.

## Pollution Prevention and Waste Management

For information on pollution prevention and waste management, refer to Section 3.

## REQUIRED REAGENTS AND STANDARDS

Description	Quantity Required		Cat. No.
	per test	Unit	
Potassium 1 Reagent Powder Pillows .....	1 pillow .....	25/pkg .....	14321-98
Water, deionized .....	10 mL .....	4 liters .....	272-56

## REQUIRED EQUIPMENT AND SUPPLIES

DR/4000, 1-inch cell, Adapter .....	1 .....	each .....	48190-00
Clippers, for opening powder pillows .....	1 .....	each .....	968-00
Sample Cells, matched pair, 1-inch, glass, with stopper .....	2 .....	pair .....	26126-02

## OPTIONAL REAGENTS AND STANDARDS

Nickel Standard Solution, 1000-mg/L .....	100 mL .....	14176-42
Nickel Standard Solution, 1000-mg/L .....	500 mL .....	14176-49

## OPTIONAL EQUIPMENT AND SUPPLIES

Filter Paper, folded, 12.5-cm .....	100/pkg .....	1894-57
Funnel, analytical poly, 65-mm .....	each .....	1083-67



FOR TECHNICAL ASSISTANCE, PRICE INFORMATION AND ORDERING:  
In the U.S.A. – Call toll-free 800-227-4224  
Outside the U.S.A. – Contact the HACH office or distributor serving you.  
On the Worldwide Web – [www.hach.com](http://www.hach.com); E-mail – [techhelp@hach.com](mailto:techhelp@hach.com)

HACH COMPANY  
WORLD HEADQUARTERS  
Telephone: (970) 669-3050  
FAX: (970) 669-2932