



Elexy™

Redefining possibilities



Volume



Viscosity



| | |
|-------------------------|---|
| Device Type | Reusable electromechanical autoinjector |
| Primary Container | 1.0 mL and 2.25 mL pre-filled syringe |
| Usage | Single dose Reusable power unit |
| Dosage Type | Empty-all Full contents |
| Route of Administration | Subcutaneous |
| Activation Method | Needle cover |
| Mode of connectivity | Cellular (LTE) |
| Injection Feedback | Audio, continuous visual (drug inspection window and LED light), and tactile feedback |

Handling Instructions



1
Insert PFS drug cassette into power unit, twist until dots align.



2
After top light turns on, uncap.



3
Press down on needle cover and wait until top light fades, indicating injection completion.



4
Remove drug cassette and dispose into sharp container.

Key Features

- ✓ One interface that is compatible to both pre-filled syringes and cartridges.
- ✓ Features a reusable and rechargeable power unit.
- ✓ Adaptable device for flexible drug volumes, drug viscosities, and injection speeds.
- ✓ Integrated cellular connectivity for seamless injection experience.
- ✓ Multiple injection feedback, including visual, audio, and tactile.



The SHL Medical logo and Elaxy are registered and/or non-registered trademarks of SHL Medical AG or its affiliates, in the United States, Europe, Taiwan, and other jurisdictions. Final product specifications are subject to change based on customer requirements and SHL's offerings. Contact info@shl-medical.com for more information.

www.shl-medical.com

Elexy™

Redefining possibilities



Volume 

Viscosity 



| | |
|-------------------------|---|
| Device Type | Reusable electromechanical autoinjector |
| Primary Container | 3.0 mL and 5.0 mL cartridge with Needle Isolation Technology(NIT®) |
| Usage | Single dose Reusable power unit |
| Dosage Type | Empty-all Full contents |
| Route of Administration | Subcutaneous |
| Activation Method | Needle cover |
| Mode of connectivity | Cellular (LTE) |
| Injection Feedback | Audio, continuous visual (drug inspection window and LED light), and tactile feedback |

Handling Instructions



1 Insert cartridge drug cassette into power unit, twist until dots align.



2 After top light turns on, unscrew cap.



3 Press down on needle cover and wait until top light fades, indicating injection completion.



+ Remove drug cassette and dispose into sharp container.

Key Features

- ✓ One interface that is compatible to both pre-filled syringes and cartridges.
- ✓ Features a reusable and rechargeable power unit.
- ✓ Adaptable device for flexible drug volumes, drug viscosities, and injection speeds.
- ✓ Integrated cellular connectivity for seamless injection experience.
- ✓ Multiple injection feedback, including visual, audio, and tactile.



The SHL Medical logo and Elaxy are registered and/or non-registered trademarks of SHL Medical AG or its affiliates, in the United States, Europe, Taiwan, and other jurisdictions. Final product specifications are subject to change based on customer requirements and SHL's offerings. Contact info@shl-medical.com for more information.

www.shl-medical.com