

NEXT GENERATION SDD FOR BEST EDXRF PERFORMANCE

# HighSense™ XP Detector

## Superior EDXRF Performance

---

Dr. Adrian Fiege  
Product Manager XRF  
10 June 2021

# HighSense™ XP Detector

## Superior EDXRF Performance

---

01

What is new?  
Next Generation SDD Technology

02

What are the benefits?  
Application Examples and QUANT Solutions

03

Why is Software important?  
Inside SPECTRA.ELEMENTS

04

What else is important?  
Benefits of the S2 PUMA Series 2

05

Summary

06


Q&A Session

# HighSense™ XP Detector

## Now available for the S2 PUMA Series 2

### State-of-the-Art Hardware & Next Generation Software

- New premium detector for all elements (C to Am)
- Robust, high transmission Graphene window (non-toxic)
- Bruker's detector chip technology
- Further enhanced cooling (Peltier) performance
- New SPECTRA.ELEMENTS with Dynamic Detector Profiling




# HighSense™ XP Detector

## Now available for the S2 PUMA Series 2

### Ultra-high transmission

- The 0.9 μm Graphene Window offers ultra high transmission
- Sensitivity improved by typically 30% compared to the HighSense LE
- A Si-grid is not required, minimizing diffraction peaks and improving the performance for Si



HighSense XP: Best performance for all element!

# HighSense™ XP Detector


## Now available for the S2 PUMA Series 2

	HighSense XP	HighSense LE	HighSense
Detector Type	SDD	SDD	SDD
Element Range	C – Am	(C) F – Am	(Na) Mg - Am
Window	Graphene (0.9 µm)	Polymere window	Be (8 µm)
Light element sensitivity ( $\leq$ Si)	Excellent	Excellent	Good
Medium element sensitivity ( $>$ Si)	Excellent	Very Good	Very Good
Heavier element sensitivity ( $>$ Fe)	Excellent	Very Good	Excellent
Dynamic Detector Profiling	Yes	No	No
Best Energy Resolution at Mn kA	<132 eV*	~135 eV	~135 eV
Max. input count rate	2.400.000 cps	1.500.000 cps	1.500.000 cps

\*At input count rate of max. 4000.000 cps

# HighSense™ XP Detector

## Now available for the S2 PUMA Series 2



Optimal excitation of the sample is ensured by:

- High power 50 Watt X-ray tube
- Up to 2 mA and 50 kV (30 kV)
- 10-position primary beam filter
- Direct excitation
- Closely coupled optics


HighSense turns power into performance

# HighSense XP

## Optimal Performance even for Trace Elements

### Lower Limit of Detection (LLD)

- Determined for optimal conditions
- $LLD = 3 \times \text{Standard Deviation}$



# GEO-QUANT Basic S2 PUMA Series 2 with HighSense XP

---

## Versatile Solution for Geo-Materials

- Powerful analytical solution developed for the analysis of major and minor oxides in geological materials.
- Based on 20 certified reference materials (CRM), incl. cement, clay, feldspar, limestone, dolomite, gypsum, bauxite, refractory, and more
- 4 dedicated QC samples
- 14 key elements covered






# GEO-QUANT Basic S2 PUMA Series 2 with HighSense XP

## Ready for various applications

- Mining operations
- Industrial minerals
- Raw materials for cement and building materials
- Refractories
- Ceramics and glass
- Geochemical research
- Archaeology
- Environmental studies

Compound	Number of CRMs	Max. conc. (wt%)
Na <sub>2</sub> O	16	11
MgO	18	40
Al <sub>2</sub> O <sub>3</sub>	19	85
SiO <sub>2</sub>	19	100
P <sub>2</sub> O <sub>5</sub>	16	7.5
SO <sub>3</sub>	13	58
K <sub>2</sub> O	18	12
CaO	18	100
TiO <sub>2</sub>	19	4
Cr <sub>2</sub> O <sub>3</sub>	14	0.1
MnO	16	0.9
Fe <sub>2</sub> O <sub>3</sub>	20	100
ZnO	11	0.2
SrO	13	0.3

# GEO-QUANT Basic S2 PUMA Series 2 with HighSense XP




Gypsum – concentrations in wt%

# GEO-QUANT Basic S2 PUMA Series 2 with HighSense XP

## Detector Comparison

HighSense XP vs. HighSense

- QC samples tested:
  - Gypsum
  - Iron Ore Sinter
  - Limestone
  - Cement
- HighSense XP: Advantages for low levels of Na, Mg, and Al.



# CEMENT-QUANT

## S2 PUMA Series 2 with HighSense XP

---

- Provides a workflow to meet good laboratory practice (GLP)
- Includes a quality check procedure for instrument stability and a repeatability check for the sample preparation
- Contains a set of predefined applications for the XRF analysis of process materials like raw mix, clinker, and cement



CEMENT-QUANT: Makes you ready for ASTM C114 and DIN EN ISO 29581-2

# CEMENT-QUANT Basic S2 PUMA Series 2 with HighSense XP

CEMENT-QUANT comes with a set of international certified reference materials (CRM)

- Allows to create [traceable](#) and [accurate](#) calibrations in three steps:
- Step 1 – Preparation of the CRMs as fused beads
- Step 2 – Measurement of the beads with CEMENT-QUANT
- Step 3 – Calibration - [Ready to analyze! Traceable to NIST!](#)



CEMENT-QUANT: Makes you ready for ASTM C114 and DIN EN ISO 29581-2

# CEMENT-QUANT Basic

## S2 PUMA Series 2 with HighSense XP

**Excellent repeatability!**

**Also for light elements like Na and Mg – thanks to HighSense XP**

(wt%)	SiO <sub>2</sub>	Al <sub>2</sub> O <sub>3</sub>	Fe <sub>2</sub> O <sub>3</sub>	CaO	MgO	SO <sub>3</sub>	Na <sub>2</sub> O	K <sub>2</sub> O	TiO <sub>2</sub>	P <sub>2</sub> O <sub>5</sub>	ZnO	Mn <sub>2</sub> O <sub>3</sub>	Cr <sub>2</sub> O <sub>3</sub>	SrO
Rep-1	31.31	9.37	0.965	51.51	4.53	3.36	0.171	0.68	0.69	0.027	0.014	0.239	0.005	0.080
Rep-2	31.31	9.38	0.962	51.52	4.55	3.36	0.173	0.68	0.69	0.029	0.014	0.238	0.005	0.081
Rep-3	31.28	9.38	0.963	51.51	4.55	3.36	0.156	0.67	0.69	0.028	0.014	0.238	0.006	0.080
Rep-4	31.31	9.39	0.964	51.52	4.53	3.37	0.150	0.67	0.69	0.031	0.014	0.238	0.005	0.080
Rep-5	31.31	9.40	0.964	51.52	4.54	3.36	0.156	0.67	0.69	0.027	0.014	0.239	0.005	0.080
...	...	...	...	...	...	...	...	...	...	...	...	...	...	...
Rep-20	31.27	9.37	0.964	51.52	4.54	3.37	0.157	0.68	0.69	0.032	0.014	0.238	0.006	0.080
Mean	31.29	9.37	0.962	51.51	4.53	3.36	0.155	0.68	0.69	0.028	0.014	0.238	0.005	0.080
Std.Dev	0.02	0.02	0.001	0.02	0.02	<0.01	0.018	<0.01	<0.01	0.002	<0.001	0.001	0.001	<0.001
Rel.Std.Dev	0.07%	0.17%	0.15%	0.04%	0.41%	0.11%	11.77%	0.30%	0.24%	7.55%	0.71%	0.38%	10.42%	0.39%

# SMART-QUANT FP

## Best standardless analysis with EDXRF

---

SMART-QUANT FP is set up to work in full Fundamental Parameter (FP) mode – this means no calibrations necessary! Excellent for raw material testing and whenever special samples outside the analytical routing need to be measured.


- Fluorine to Americium
- ppm to 100%
- Air, Helium, Vacuum
- 30 and 50 kV



SMART-QUANT: Push-button solution for quick and reliable analysis of unknown samples

# SMART-QUANT FP

## Powered by HighSense XP




Nickel Alloy – Waspaloy, concentrations in wt%



# SMART-QUANT FP


## Powered by HighSense XP



High Al Cement, concentrations in wt%

# The Benefits of Detector Profiling

- New, fully integrated feature in SPECTRA.ELEMENTS
- Flexible: Make your selection each analytical range
  - **High Resolution** for optimal peak separation of neighboring elements
  - **High-Count-Rate** to boost the throughput or decrease LLD




Example: Achieve 4 x higher net intensity for **Cement** samples with the High-Count-Rate profile

**SPECTRA.ELEMENTS**

Easy. Fast. Smart. Powerful.

# The Benefits of Detector Profiling

- New, fully integrated feature in SPECTRA.ELEMENTS
- Flexible: Make your selection each analytical range
  - **High Resolution** for optimal peak separation of neighboring elements
  - **High-Count-Rate** to boost the throughput or decrease LLD



Name	Voltage (kV)	Filter	Detector Profile
[-] Range 1	40	Al (500μ)	Standard
Cu KA1/test	40	Al (500μ)	Standard
[-] Range 1/40kV	40	Al (500μ)	High resolution
Zn KA1/test	40	Al (500μ)	High resolution
[-] Range 1/40kV/Al (5...	40	Al (500μ)	High count rate
Ni KA1/test	40	Al (500μ)	High count rate

Select the profile for each range when you setup your method.

## SPECTRA.ELEMENTS


Easy. Fast. Smart. Powerful.

# The Benefits of Detector Profiling

## Precision Test

- Cement QC sample
- 3 ranges, 190 or 600 s total counting time, 25 repetitions

The [High-Count-Rate](#) setting allows to reduce the counting time factor of ~3




# The Benefits of Detector Profiling

## ISO 29581 / EN 196-2:

A Cement norm which defines two categories of precision (repeatability):  
“Normal” and “Expert”

- The S2 PUMA Series 2 with HighSense XP meets standard deviation limits easily.



## What requirements for a modern spectrometer software?

- Intuitive user interface
- Full functionality
- Different user levels
- Short processing times
- High stability
- Flexible data management

Modern XRF software offers better analytical performance and **saves time!**

(training, operation, evaluation)



**SPECTRA.ELEMENTS:** Next Generation XRF Software Platform

# SPECTRA.ELEMENTS

## Next Generation Spectrometer Software

Faster. Smarter. Easier.

Faster operation

- Immediate data processing
- Smooth navigation
- Rapid access to database

Solution				Preparation		
Solution		SMART-Oxides (201)		Values	Preparation = Fused bead (ig...	
Inputs				Basis	Original sample	
	Z	Formula	Concentration	Evaluation ...	Quantified by	Range Name
▶	13	Al <sub>2</sub> O <sub>3</sub>	1,39 %	Stdless	None	50 kV, Range 1
	14	SiO <sub>2</sub>	35,88 %	Stdless	None	50 kV, Range 1
	16	SO <sub>3</sub>	9,16 %	Stdless	None	50 kV, Range 1
	17	Cl	2,54 %	Stdless	None	50 kV, Range 1
	19	K <sub>2</sub> O	0,08 %	Stdless	None	50 kV, Range 1
	21	Sc <sub>2</sub> O <sub>3</sub>	0,03 %	Stdless	None	50 kV, Range 1
	22	TiO <sub>2</sub>	0,11 %	Stdless	None	50 kV, Range 1
	25	MnO	0,06 %	Stdless	None	50 kV, Range 1
	26	Fe <sub>2</sub> O <sub>3</sub>	0,12 %	Stdless	None	50 kV, Range 1
	27	CoO	0,00 %	Stdless	None	50 kV, Range 1
	38	SrO	0,01 %	Stdless	None	50 kV, Range 1
	39	Y <sub>2</sub> O <sub>3</sub>	0,00 %	Stdless	None	50 kV, Range 1
	41	Nb <sub>2</sub> O <sub>5</sub>	0,00 %	Stdless	None	50 kV, Range 1
	45	Rh <sub>2</sub> O <sub>3</sub>	0,00 %	Stdless	K	50 kV, Range 1
	82	PbO	0,00 %	Stdless	None	50 kV, Range 1
	20	CaO	50,61 %	Stdless	None	
*						

Evaluate your measurements in no time!


# SPECTRA.ELEMENTS

## Next Generation Spectrometer Software

Faster. Smarter. Easier.

Smarter features

- Dedicated User Levels
- Basic & Advanced Mode
- Detector Profiling
- Aut-O-Matrix
- Loss-Eliminated Alphas,
- Pressure compensation
- Blanks, duplicates, QC
- Post-Processing
- Automation ready





# SPECTRA.ELEMENTS

## Next Generation Spectrometer Software

Faster. Smarter. Easier.

Easier to use

- Quick learning with new User Interface
- Hints provide help when needed
- WIAZARD: the proven tree-structure guide you through the calibration process
- LOADER: Intuitive interface for routine operation
- RESULT MANAGER: Access all your data quickly and get extended reporting



The LOADER: load / unload samples; adjust priorities; start / stop your measurements; view your results.

# S2 PUMA Series 2

## Taylor-made Solution for all Applications

### Benefits of the S2 PUMA

- Sample Handling for all needs:
  - Single
  - XY-Autochanger
  - XY-Automation
  - Carousel
  - Mapping-Stage



# S2 PUMA Series 2

## Taylor-made Solution for all Applications

---

### Benefits of the S2 PUMA

- All atmospheres
  - Air, Vacuum, Helium, Nitrogen.
- All samples
  - Liquid, Powder, Pressed Pellet, Solid, Bulk.




# S2 PUMA Series 2

## Taylor-made Solution for all Applications

### Benefits of the S2 PUMA

- Easy-of-Use
  - SPECTRA.ELEMENTS
  - TouchControl™
- High uptime and short maintenance
  - SampleCare™
  - Smart maintenance features
  - Self-diagnostic functions
  - Remote service & support




# Summary

## S2 PUMA Series 2 with HighSense XP



- HighSense XP combines ultra-high sensitivity with robustness
- Detector Profiling allows to reduce counting times or decrease LLDs significantly
- SPECTRA.ELEMENTS is now even faster and comes with additional features to make your life easy (e.g. Basic and Advanced mode)
- Ready to perform from day one with out-of-the-box solutions
  - SMART-QUANT, GEO-QUANT, CEMENT-QUANT






# Q&A Session

---

Dr. Adrian Fiege  
[Adrian.Fiege@bruker.com](mailto:Adrian.Fiege@bruker.com)





Innovation with Integrity