
HardyCHROM[™] NG-Test[®] CARBA 5



Where **EASY-TO-READ** chromogenic culture media meets RAPID, CLEAR RESULTS to combat ANTIMICROBIAL RESISTANCE.

Used TOGETHER, you can fight back against the global threat of CARBAPENEM-RESISTANT **ENTEROBACTERALES (CRE).**



Often called the "nightmare bacteria" CRE are a family of germs that are difficult to treat because they are resistant to nearly all antibiotics, including carbapenems, some of our most powerful drugs.

Rapid spread of CRE in parts of the United States and other countries in both humans and animals points to a need to more aggressively detect and control CRE spread.



Estimated cases in hospitalized patients



1,100 Estimated deaths




\$130M Estimated attributable healthcare costs in

(Source: 2019 AR Threats Report)

CARBAPENEM- RESISTANT Enterobacterales


With these two tools, you are in the fight against the worldwide threat of CRE.



Using HardyCHROM™ CRE + CARBA 5, saves time, money and resources, but more importantly, accelerates detection for faster and more accurate treatment.

NG-Test® CARBA 5

NG-Test® CARBA 5 is the only rapid, multiplex, phenotypic test capable of detecting KPC, OXA-48-like, VIM, IMP, and NDM carbapenemases produced by Enterobacterales and Pseudomonas aeruginosa.




CLEAR RESULTS HELP TO GUIDE THERAPY AS WE COMBAT THE WORLD'S DEADLIEST FORCES






User Friendly - 5 simple steps!



A "LINE" means an ENZYME!

Example of interpretation



Harc CHROM

FOR THE SCREENING OF ESCHERICHIA COLI. KLEBSIELLA AEROGENES, KLEBSIELLA OXYTOCA, KLEBSIELLA PNEUMONIAE, ENTEROBACTER CLOACAE COMPLEX, AND SERRATIA MARCESCENS.

HardvCHROM™ CRE is a selective and differential chromogenic agar medium intended for the qualitative and presumptive detection from stool specimens of Escherichia coli that are non-susceptible to carbapenems as pink colonies and KES (Klebsiella aerogenes, Klebsiella oxytoca, Klebsiella pneumoniae, Enterobacter cloacae complex, and Serratia marcescens) that are non-susceptible to carbapenems as blue colonies.


HardyCHROM™ CRE is intended as an aid in the detection, identification of colonization and control of these bacteria in a healthcare setting. HardyCHROM™ CRE is not intended to diagnose infection or guide therapy. Results can be interpreted after incubation for 18-24 hours. Subculture to non-selective medium is required for confirming identification, antimicrobial susceptibility testing and epidemiological typing.



15x100mm, 10/pk Cat. no. G323


Colonies that are not pink to magenta, blue, or blue with pink halos are negative; no carbapenem non-susceptible Escherchia coli or KES detected.

EASY TO READ RESULTS



Non-susceptible *Escherichia coli* produces colonies that are rose to magenta in color with darker pink centers.

Klebsiella, Enterobacter, and Serratia spp. produce large, dark blue colonies (with or without pink halos) are presumptive positive for carbapenem non-susceptible KES (Klebsiella aerogenes, Klebsiella oxytoca, Klebsiella pneumoniae, Enterobacter cloacae complex, and Serratia marcescens).











FM 572526

Hardy Diagnostics adheres to cGMP, is licensed by the FDA as a medical device manufacturer, and its quality management system is ISO 13485 certified.


Headquarters

1430 West McCoy Lane Santa Maria, CA 93455 800.266.2222 Sales@HardyDiagnostics.com HardyDiagnostics.com

Distribution Centers

Santa Maria, California Olympia, Washington Salt Lake City, Utah Phoenix, Arizona Dallas, Texas Springboro, Ohio Lake City, Florida Albany, New York Raleigh, North Carolina

Copyright © 2022 Hardy Diagnostics



FREE SAMPLES AVAILABLE!

For more information, including quotes, free product samples, or demonstrations, please contact us:

800.266.2222 HardyDiagnostics.com/carbapenemase-detection