

STREP B
CARROT BROTH™
ONE-STEP



 **HARDY**
DIAGNOSTICS
A Culture of Service™

Simple detection of GBS in **one step!**

Every life deserves an opportunity to **grow** and **thrive**.

Better tests enable you to provide **just that**.

**Better testing.
Better processes.
Better outcomes.**

Our exclusive Carrot Broth™ One-Step system allows for **faster results** for the **cultivation** and **identification** of beta-hemolytic strains of Group B Streptococcus in pregnant women.





Compared to conventional methods, it's not even close

When compared to methods such as PCR or LIM broth, Carrot Broth significantly reduces the work load, the cost, and time necessary to accurately identify GBS colonization in pregnant women.

Strep B Carrot Broth™ One-Step features greater sensitivity and decreased turn-around-time, which lowers overall costs.

A clear, chromogenic reaction as early as 16 hours after testing tells you everything you need to know: **Positive.**

Positive results do not require additional work ups.

The advantages are clear

- Easy one-step process! Place the patient specimen directly into the media.
- Development of any orange to red color indicates a positive result.
- For positives, there is no need for further subculturing to confirm.
- Significantly lower costs when compared to PCR methods.
- **Increased accuracy when compared with LIM Broth. Clinical testing demonstrated a sensitivity of 98.8% and a specificity of 98.2% for the recovery of GBS from patient specimens when compared to LIM Broth.***
- **Strep B Carrot Broth™ One-Step** was found to be more sensitive than PCR!***
- Negative cultures from **Strep B Carrot Broth™ One-Step** should be subcultured to **GBS Detect™** (Cat. no. A300) to detect weakly or non-hemolytic strains of GBS.





- **Strep B Carrot Broth™ One-Step** is available in multiple formats that work on many automated inoculating platforms.
- In house studies showed a sensitivity of 99% at 10² CFU/ml, and specificity of 100%. Data available upon request.
- *Ledeboer, N. et al. **“Evaluation of New Strep B Carrot Broth™ One-Step in the Detection of Group B Streptococcus: A Multi-Center Study”** poster at ASM 2017, New Orleans, LA. During the study, there were 9 beta hemolytic GBS recovered by color reaction in Carrot Broth that were not recovered by LIM broth. There were also 4 non-hemolytic GBS recovered upon subculture of Carrot Broth to blood agar that were not recovered by LIM broth.
- **Schreckenberger, P. **“Comparison of Real-time PCR with standard LIM Broth and Carrot Broth for the Detection of Group B Strep”** poster at ASM 2006, Orlando, FL.

To further explore the research listed above, please visit:
hardydiagnostics.com/carrot-broth/

Strep B Carrot Broth™ One-Step	
13x100mm glass tube,20/pk	Z40
Strep B Carrot Broth™ One-Step	
12x80mm plastic tube,100/pk	Z44BX
Strep B Carrot Broth™ One-Step	
16x100mm glass tube,100/pk	Z46BX
StrepPRO™ Streptococcal Grouping Kit	
A latex agglutination test for the rapid identification of beta-hemolytic streps includes groups A, B, C, D, F, and G. 60 tests	PL030HD


GBS Detect™

for Weakly or Non-Hemolytic Strains

A perfect companion product to Carrot Broth™ which does not detect non-hemolytic strains.

GBS Detect™ is recommended for the isolation and detection of weakly or non-hemolytic Group B Streptococci.

GBS Detect™ induces beta-hemolysis of weakly or non-hemolytic GBS upon subculture from enrichment broth procedures, such as LIM Broth or Strep B Carrot Broth™ One-Step.



GBS Detect™

15x100mm plate, 10/pk

A300

GBS Detect™

15x100mm plate, 100/bx

A300BX

GBS Detect™

15x100mm plate with reduced stacking ring for automated inoculating devices, 10/pk

GA300



Control Organisms

***Streptococcus agalactiae* derived from ATCC® 12386™**

KWIK-STIK™ 6/pk



0439K

***Streptococcus pyogenes*, Group A derived from ATCC® 19615™**

KWIK-STIK™ 6/pk

0385K

Strep B Carrot Broth™ One-Step method





HARDY DIAGNOSTICS

A Culture of Service™

4
MORE THAN



Years of
Bringing You
A Culture of Service



100%
Employee Owned

ESTABLISHED
1980

MADE IN
USA



Hardy Diagnostics
donates 1% of net
profits to charity.



FM 572526

Hardy Diagnostics has a Quality Management System that is certified to ISO 13485 and is a FDA licensed medical device manufacturer.

Headquarters

1430 West McCoy Lane
Santa Maria, CA 93455
800.266.2222
Sales@HardyDiagnostics.com
HardyDiagnostics.com

Distribution Centers

Santa Maria, California
Olympia, Washington
Salt Lake City, Utah
Phoenix, Arizona
Dallas, Texas
Des Moines, Iowa
Springboro, Ohio
Lake City, Florida
Albany, New York
Raleigh, North Carolina



Contract Holder

Contract # 47QSWA22D0077