

# START HERE GET TO KNOW THE BUTTONS & ATTACHMENTS

Charge your unit for 3 hours prior to first use.

Please carefully review all Warnings, Precautions, and Contraindications in User Guide prior to first use.

**EXFOLIATION ATTACHMENT:**

**INTENSITY BUTTON**  
Short press to change suction level.

**POWER BUTTON**  
Short press to pause and resume.

**DEPUFFI™ ATTACHMENT:**

**INTENSITY BUTTON**  
Short press to change temperature intensity.

**POWER BUTTON**  
Short press to switch between INSTACHILL and INSTAHEAT settings.

## EXFOLIATION BENEFITS

- ENHANCED GLOW**  
Radiant, lit-from-within glow.
- SUPERCHARGED EXFOLIATION**  
Combines AHA/BHA and physical exfoliation for brighter, more radiant skin.
- POWERED DETOX**  
Clearer, less visible pores.
- BOOSTED HYDRATION**  
Leaves skin deeply hydrated for a smoother feel.

## DEPUFFI™ BENEFITS

- INSTACHILL**  
**ADVANCED DEPUFFING**  
Reduces appearance of dark circles and leaves under eyes refreshed and visibly less puffy.
- INSTACHILL**  
**SOFTENING PREP**  
Primes skin for extraction by loosening surface impurities.
- INSTACHILL**  
**RAPID SOOTHING**  
Cooling relief for instantly calmer soothed skin.
- INSTAHEAT**  
**ACTIVATED RADIANCE**  
Massage-driven circulation for a healthy, radiant look.
- INSTAHEAT**  
**FACIAL DE-STRESS**  
For muscle tension relief.

# CLINICALLY TESTED SKINCARE

Dermatologist-Tested Korean skincare formulas.  
Developed for optimal results when using FacialPro™ Glow device.



## DERM DETOX

**AHA + BHA Exfoliating Gel**  
Gently softens and dissolves excess dead skin cells and sebum sitting on surface of skin.

APPLY



## HYDRO INFUSE

**BHA Hydrator**  
Hyaluronic acid infuses skin with moisture.

INFUSE



**Pro Tip:** Don't forget to use both Wide and T-Zone tips when using your Hydro Boost for best at-home results.

TO REPLENISH SKINCARE, VISIT [SHARKBEAUTY.COM](https://www.sharkbeauty.com)

Only use this skincare with FacialPro™ Glow system.



SCAN TO REPLENISH

# Shark FACIALPRO™ GLOW

## + DePuffi

### QUICK START GUIDE



# DAILY DEPUFFING ROUTINE


Depuff to sculpt and boost radiance daily, in just 4 minutes.

## GET STARTED

- Clean your face to remove any dirt and oil.
- Connect Depuffi™ Attachment.
- Apply your favorite serum, face oil or moisturizer for best glide.


## INSTACHILL / 4 MIN. SESSION

- Turn the device on.
- Device starts at temperature level 3. Start there and dial down to your preferred temperature setting.
- Glide over face, neck, jawline, and shoulders to reduce puffiness for a sculpted look and refreshed skin.



## INSTAHEAT / 4 MIN. SESSION

- Turn the device on.
- Short press to switch into InstaHeat mode.
- Device starts at temperature level 3. Start there and dial down to your preferred temperature setting.
- Massage and glide over face, neck and jawline to boost circulation, encourage lymphatic drainage and boost skin's radiance.



## DEFINE AND DEPUFF

### INSTACHILL

- Move in slow upward and/or outward strokes to reduce puffiness for a defined jawline and cheekbones that look sculpted.
- Use under eyes to reduce the appearance of dark circles for an instantly refreshed and awakened look.



## BOOST RADIANCE

### INSTAHEAT

- Massage in slow upward and/or outwards strokes to boost circulation for more radiant skin.
- Use in downward motion on the neck for lymphatic drainage.



# WEEKLY PRO GLOW ROUTINE

Smooth, brighten, and hydrate with a weekly skin refresh and glow in just 9 minutes.



**DEPUFFI RECOMMENDATION:** Prior to extraction, prep skin with **INSTAHEAT** to loosen surface impurities.

## 1 APPLY DERM DETOX / 3 MIN. SESSION

### AHA + BHA gently dissolves impurities

- Apply a thin layer of **Derm Detox** all over face, avoiding eyes and mouth.
- Leave on for 3 minutes and proceed to next step.



## 2 EXFOLIATE & EXTRACT / 3 MIN. SESSION

### Clears pores and lifts impurities.

- Fill tank with cool water to the "MAX" fill line.
- Attach **Wide Tip** and turn device on.
- Adjust intensity and glide over skin as directed.
- Pause device and attach **T-Zone Tip** for targeted care. Device will resume at time you left off.
- Device pauses after 3 minutes—move to next step.
- Rinse off any remaining Derm Detox from face with water.



## 3 HYDRO INFUSE / 3 MIN. SESSION

### Locks in moisture and refreshes skin.

- Fill tank with **Hydro Infuse** to "MAX" fill line.
- Select your **tip** and intensity level as Step 2 and glide over skin as directed.
- Once you have completed your Weekly Pro Glow Routine, apply moisturizer and/or SPF.
- If using Depuffi™ on InstaChill as recommended below, apply moisturizer beforehand and SPF after completing the step.



**DEPUFFI RECOMMENDATION:** After extraction, refresh skin by applying your moisturizer, using **INSTACHILL** to soothe for instantly calmer skin. Complete routine by applying SPF if it's daytime.

# CHOOSING A SET OF TIPS

SELECT EITHER THE NORMAL TIPS OR THE GENTLE TIPS BASED ON YOUR PREFERENCE.

NORMAL



stronger feeling of suction  
FOR MOST SKIN TYPES

CHOOSE  
A SET



GENTLE



more mild feeling of suction  
FOR SENSITIVE SKIN TYPES

All tips were designed with performance in mind.

Dermatologist-tested to deliver clean pores and hydrated, glowing skin.

DO NOT HOLD THE TIPS STILL ON THE SKIN. KEEP THEM MOVING.




## HOW TO GLIDE WITH WIDE TIP

For use on cheeks, jawline, and upper forehead.

We recommend you use the Wide Tip for 2 out of the 3 min from both Step 2 and 3.

- Place the tip on your cheek, jawline, or upper forehead and press gently to create a seal.
- You'll feel suction and hear a soft "kiss" with each stroke.
- Move in short, smooth strokes about 1-2 inches long.
- Use your free hand to gently stretch skin for a smooth glide.




## HOW TO GLIDE WITH T-ZONE TIP

For use on T-zone: nose, chin, and around brow bones.

We recommend you use the T-Zone tip for 1 out of the 3 min from both Step 2 and 3.

- Place the more precise tip on your nose, chin, and between eyebrows and press gently to create a seal.
- You'll feel suction and hear a soft "kiss" with each stroke.
- Move in short, smooth strokes about 1-2 inches long.
- Use your free hand to gently stretch skin for a smooth glide.



The combination of both **Wide** and **T-Zone Tip** is what will ensure optimal results across all areas of the face.