currentbody.com

CurrentBody.com Ltd. Q17, Crossley Road, Stockport, Greater Manchester,
SK4 5BB, United Kingdom.

CurrentBody BV, Naritaweg 2, 1437EL, Rozenburg, Netherlands

Manual Version: V1.0
Issue Date: 2022-03-08

CURRENTBODY | skin™

LED TEETH WHITENING KIT

USER MANUAL

Includes important safety information

Read all instructions prior to use

PLEASE KEEP THESE INSTRUCTIONS SAFE

Welcome to CurrentBody Skin™ LED Teeth Whitening Kit

Established in 2009, CurrentBody is the world's only online retailer dedicated to home-use beauty devices. As the first company to recognise devices as one of the fastest growing trends in the world of beauty, we have unrivalled expertise in selecting at-home devices that are science-backed and clinically proven. We offer our customers the ultimate choice of devices that treat a wide range of concerns, including anti-ageing, cleansing and hair removal.

Underpinning our range of innovative beauty technology is expert education to help our customers choose the right device for them and excellent customer service.

Our highly trained beauty device experts are on hand to offer impartial advice and answer your questions in person, as well as how-to videos and top tips available throughout our website and social media channels.

WHAT IS IN THE BOX

CurrentBody Skin™ LED Teeth Whitening device
User Manual
USB charging cord
Storage bag
Whitening Gel (3 x 5ml syringes)
Teeth shade guide

CONTRAINDICATIONS

Do not use this device if your skin or gums are sensitive to light or you are currently using medication that may cause skin sensitivity to light. Doing so may result in a photosensitive reaction.

If you are unsure whether your medication causes sensitivity to light, confirm with your physician prior to using this device.

Discontinue use immediately if you experience any sensitivity to using this device.

For further information on medication that can cause light sensitivity, please visit
<https://www.currentbody.com/pages/led-contraindications>

WARNINGS

Do not immerse this device. Keep water away from the main body of the device and charging base.

Stop using device immediately if you think it is defective.

Do not attempt to take the device apart. Any modification or attempted service will void the warranty.

Do not store the device in a sunny location with high temperatures as this may cause damage to plastic casing.

Discontinue use and contact Customer Service if device overheats or becomes uncomfortably hot on your skin or gums.

Always unplug the device from the outlet after charging is complete.

Use only the supplied Universal USB charging cord. Use of an incorrect power cord could damage the device.

Do not use if your skin or gums are swollen from a previous treatment.

Do not use the device to treat skin or gums that are painful to touch, raised (higher than surrounding skin), crusting, or bleeding.

This may be a sign of a more serious condition and should be checked by your doctor.

This device contains a Lithium Ion Battery. Before disposal, please review your local laws and requirements surrounding Lithium Ion Battery disposal. The preferred method of disposal is through recycling.

IMPORTANT SAFETY INFORMATION

This device is designed for use on adults over 18 years of age and not recommended for use on infants or children.

Do not shine the device lights directly into the eyes. Doing so may cause temporary vision discomfort. This device also contains blue light which can be harmful to eyes - ensure you place the device in the mouth before switching on.

Contact your doctor if you have any medical questions regarding use of this device.

Consult your physician prior to using this device if you:

1. Are pregnant or suspect you may be pregnant.
2. Have suspicious cancerous lesions.

Do not use or charge the device with a power adapter that is damaged or faulty. Do not harshly bend or abuse the USB charging cord.

Do not immerse the device in water and keep water away from the head of the device.

This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.

Children should be supervised to ensure that they do not play with the appliance.

CurrentBody Skin™ LED Teeth Whitening Kit

CurrentBody Skin™ LED Teeth Whitening Kit uses light emitting diodes (LEDs) in the blue, red and infrared spectrums which are clinically proven to increase circulation and enhance gum health.

CHARGING INSTRUCTIONS

1. There are 2 ways to charge the device:

- **From the charging base:**

Place the device on the base, ensuring the charging points on either side are in contact with the base. Plug the USB charging cord into the charging base. Plug the USB charging cord into a Universal USB Power Adapter (not included) and then into a mains outlet. The power indicator light will flash when charging.

- **From the USB charging point on the device:**

Plug USB charging cord into the USB port connector on the device. Plug USB charging cord into a Universal USB Power Adapter (not included).

2. Charge your device for 2-4 hours.

3. Device is fully charged when indicator light is solid white.

LIGHT-SENSITIVITY TEST

Test yourself for light-sensitivity prior to initial treatment:

1. Power on your device.
2. Place device on of the skin on your forearm for three minutes.
3. After 3 minutes, check the area. If area has turned red, and redness lasts for more than 2 hours, your skin is light-sensitive. Please return the device to CurrentBody within 14 days for a full refund.

HOW TO USE

1. Brush your teeth.
2. Using the shade guide, match your teeth to the nearest shade colour.
3. Apply a thin layer of whitening gel to the front of the top and bottom mouth trays, using a quarter of a syringe on each tray (half a syringe in total).
4. Place silicone tray in your mouth. It is important you place the tray in your mouth before you switch the device on.
5. Press and hold the ON button until the light comes on.
6. Use the device for the 20 minute teeth whitening treatment. After 20 minutes the device will automatically switch off. You can switch the device off at any point during the treatment by pressing the ON / OFF button.

7. Rinse your mouth with water after use.
8. Please charge the device after each use. Once fully charged, unplug the device. Leaving the device on charge when not needed will shorten the life of the battery.
9. Place device into protective storage case when not in use.
10. After using for 6 consecutive days, check your teeth against the shade guide to see the improvement. The average improvement is 6 shades over 6 days. If you require further whitening, continue to use daily or use as required for top up sessions.

MAINTENANCE AND CLEANING

Maintenance:

1. This device has no serviceable parts. The battery is not accessible and not serviceable. Users should not attempt to access or service the battery as this will damage the device and void the warranty.
2. Device includes only an ON/OFF button. All settings are set to the factory setting.
3. Charge the device only with the supplied USB power charging cord and a universal USB power adapter (not included).
4. The power indicator light will begin to flash when battery levels are low. This indicates the device should be recharged.
5. Compatible Power Supply: 100V/240V

Cleaning:

1. Turn device off and unplug from any power source.
2. Separate the device from charging base.
3. Wash only the mouth tray in luke warm water. Avoid getting water on the rest of the device. Do not immerse the device in water. Keep water away from the body of the device and base.
4. Ensure the mouth tray is completely dry before storing.
5. To remove stubborn stains from the mouth tray, brush with a small amount of toothpaste and rinse.

TROUBLESHOOTING

System is not turning on:

1. Is the device fully charged? (power indicator light is on and solid)
2. Is the power indicator light illuminated?

System is not charging:

1. Is USB power cord firmly inserted into the device?(does power indicator light turn on?)
2. Try a different wall outlet to determine if there is power to the outlet.

System shuts off abruptly, prior to completion of 20 minute treatment:

1. System is designed to shut off automatically after 20 minutes.
2. System may shut off after multiple repeated uses if the device exceeds desired temperature. Allow the device to cool for several minutes and try again.

Some lights appear turned off:

1. The infrared (IR) light spectrum is not visible to the human eye. These lights will look as if they are not on, but they are operating properly. If the red, visible lights are on, then the device is operating as designed.

CUSTOMER SERVICE

Please email sales@currentbody.com if you have any questions or you need any advice.

WARRANTY

Under normal use this product is guaranteed for 2 years after its original purchase date against manufacturing defects in materials or workmanship. This warranty extends only to the original purchaser. If the product should become faulty within the warranty period, contact CurrentBody.

- CurrentBody reserves the right to repair a defective product.
- CurrentBody reserves the right to replace a defective product with the most compatible product currently available.
- This warranty does not cover damage caused by the following:
 - Accident, misuse, abuse or alteration
 - Servicing by unauthorised persons
 - Connection to incorrect power sources
 - Use with unauthorised accessories

CurrentBody shall not be responsible for any incidental, special or consequential damages resulting from the use of this product. All implied warranties, including but not limited to implied warranties of fitness and merchantability are limited in duration to two years from date of original purchase. This warranty gives you specific legal rights, and you may also have other rights that vary from country to country. Some countries do not allow the exclusion or limitation of incidental or consequential damages, or limitations on how long an implied warranty lasts, so some of the above limitations and exclusions may not apply to you. This warranty is given in lieu of all other warranties, written or oral, whether expressed by affirmation, promise, description, drawing, model or sample. Any and all warranties other than this one, whether expressed or implied, including implied warranties of merchantability and fitness for a particular purpose are hereby disclaimed. If the product is defective within the warranty period please contact CurrentBody.

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