



**Molly®**  
Modular Platform  
Technology



Volume   
Viscosity 

Device Type	Disposable autoinjector
Primary Container	1.0 mL long pre-filled syringe with staked needle
Dosage Type	Fixed   Full contents
Route of Administration	Subcutaneous
Activation Method	Needle cover
Injection Feedback	Audible click at beginning and toward end of injection, viewing window

Copyright © 2025 SHL Medical AG. All rights reserved. Version 2.6 Information accurate at time of print and is subject to change without notice. Contact [info@shl-medical.com](mailto:info@shl-medical.com) for more information.

## Handling Instructions



Remove cap.



Press and hold against injection site to activate.



Discard in sharps container.

## Key Features

- ✓ Simple 2-step operation
- ✓ Needle cover activation
- ✓ Audible, visual, and tactile feedback
- ✓ One-handed operation
- ✓ Ultra-compact design
- ✓ Permanently hidden needle
- ✓ Passive locking needle cover
- ✓ Preconfigured or customized to client specifications