



# MAGNITUDE™

When more matters



Volume



Viscosity



|                         |                                                              |
|-------------------------|--------------------------------------------------------------|
| Device Type             | Disposable autoinjector                                      |
| Primary Container       | 3.0 mL and 5.0 mL cartridges with pre-installed needle (NIT) |
| Dosage Type             | Fixed dose   Optional automated priming and titration        |
| Route of Administration | Subcutaneous                                                 |
| Activation Method       | Needle cover                                                 |
| Injection Feedback      | Continuous audible clicks during injection, viewing window   |

Copyright © 2026 SHL Medical AG. All rights reserved. Version 1.0 Information accurate at time of print and is subject to change without notice. Contact [info@shl-medical.com](mailto:info@shl-medical.com) for more information.

## Handling Instructions

\*For suspension drugs that require a priming step, please hold device with cap upright during cap removal.



1  
Unscrew cap until needle cover automatically extends\*. Then remove the cap.



2  
Press and hold against injection site to activate.



+  
Discard in sharps container.

## Key Features

- ✓ Compact design
- ✓ Simple 2-step operation
- ✓ Pre-installed Needle Isolation Technology (NIT)
- ✓ Permanently hidden needle
- ✓ Customizable needle length and gauge
- ✓ Customized to client specifications