



# EvoEndo Model LE 110 Gastroscope



We achieve  
**big things**  
through  
**small scopes.**

The EvoEndo Model LE 110 Gastroscope is the only endoscopy product FDA-cleared for both **pediatric and adult patients**, introducing a sedation-free, lower-risk, cost-effective alternative to traditional upper endoscopy.

Our ultra-slim technology unleashes the potential of sedation-free endoscopy to improve the lives of patients, families, and healthcare professionals.



## Benefits

### + Safe

- The ONLY ultra-slim single-use gastroscope cleared by the FDA specifically for Upper GI endoscopy in patients of all ages.
- Enables Sedation-Free TNE, eliminating anesthesia risk.
- Single-use design minimizes cross-contamination.

### + Capable

- 2.0 mm channel for standard pediatric accessories.
- 4-way steering with air, water, and suction functionality.

### + Flexible

- Portability makes system set-up a breeze.
- Low start-up costs make Sedation-Free TNE possible anywhere.
- Maximizes workflow efficiencies.





## EvoEndo Controller

Lightweight and portable at 6" x 8.5" and 2 lbs.




## Enhance Image Feature

Unmatched definition and clarity in a single-use disposable scope.



## Ergonomic Design

Familiar control layout fits a range of hand sizes for precision and confidence.



## Optical System

Field of View	120°
Illumination method	LED

## Insertion Portion

Bending section	210° Up, 90° Down, 100° Left, 100° Right
Maximum diameter of insertion portion	3.5 mm (0.14")
Working length	110 cm

## Connections

Air Connector	Connects to 1/4" – 3/8" supply lines
Water Connector	Connects to bottle with 1.25" top
Suction Connector	Connects to 6 mm suction device

## Channel

Average inner diameter	2.0 mm (0.078")
------------------------	-----------------



## Ready to learn more?

Scan the QR code to visit our website!

Discover clinical evidence, instructional videos, schedule a hands-on demonstration, and more!



## Patient-Centric

The EvoEndo System is designed to make endoscopy more convenient and accessible, to shorten time to diagnosis, and to eliminate the need for anesthesia.

Our ultra-slim scope accommodates small anatomy, and our Patient Experience Kit includes VR goggles for patient distraction.



EMAIL [info@evoendo.com](mailto:info@evoendo.com)

TEL 303.223.7445

FAX 303.265.9381

WEB [evoendo.com](http://evoendo.com)

 [evoendo](https://www.linkedin.com/company/evoendo)

 [evoendoinc](https://www.facebook.com/evoendoinc)

 [evoendo](https://twitter.com/evoendo)

 **evoendo**<sup>®</sup>  
Evolving Endoscopy

\*Requires Android OS mobile device with active wi-fi connection

© EvoEndo Inc.® US: Rx only. For use by trained physicians only. MKT-154 Rev 1.0