

Glossary of Terms

CPT	Current Procedural Terminology
EGD	Esophagogastroduodenoscopy. This is also known as upper gastrointestinal endoscopy.
IFU	Instructions for Use
TNE	Transnasal Endoscopy
TN-Eso	Transnasal Esophagoscopy
TN-EG	Transnasal Esophagogastroscopy
TN-EGD	Transnasal Esophagogastroduodenoscopy
VR	Virtual Reality
EvoEndo®	Medical device company developing the EvoEndo® Single-Use Endoscopy System designed for unsedated endoscopy.
EvoEndo® Comfort Kit	A specialized kit designed for patient distraction during an unsedated transnasal endoscopy.
EvoEndo® Controller	The 6"x8", 2 lbs box that powers and enables the function of the EvoEndo® Model LE Single-Use Gastroscope.
EvoEndo® Learning Series	A series of instructional videos designed to enable implementation of unsedated transnasal endoscopy at a new medical center.
EvoEndo® Model LE Single-Use Gastroscope	Sterile single-use gastroscopes with 3.5 mm outer diameter, 1.1 m length, 2 mm working channel, 4-way steering, and integrated air, water, and suction.
EvoEndo® Educator	EvoEndo's online educational resource including checklists, handouts, and videos on starting unsedated endoscopy with EvoEndo's System.