



# Evolving Sedation-Free Endoscopy

The EvoEndo® Single-Use Endoscopy System allows a full transnasal EGD exam to be performed without general anesthesia or sedation.

## Indications for Use

The EvoEndo® Model LE Gastroscope is intended for the visualization of the upper digestive tract in adults and pediatric patients, specifically for the observation, diagnosis, and endoscopic treatment of the esophagus, stomach, and duodenal bulb in patients over the age of five. The gastroscope is a sterile, single-use device and can be inserted orally or transnasally. The EvoEndo® Controller is intended for use with an EvoEndo® Endoscope for endoscopic diagnosis, treatment, and video observation.

## EvoEndo® Model LE Single-Use Gastroscope

- 3.5 mm outer diameter
- 2 mm working channel
- 1.1 m length
- HD camera
- 4-way steering
- Integrated air, water, suction

## EvoEndo® Controller

- 6" x 8.5"
- 2 lbs

## EvoEndo® Patient Experience Kit

- Distraction and comfort items to improve the patient experience



© EvoEndo, Inc. US: Rx only. For use by trained physicians only. For complete product information, warnings, and instructions, please read our instructions for use at [evoendo.com/ifu](http://evoendo.com/ifu).

# EvoEndo® Education and Training

We provide documentation, videos, and tools to train, guide, and support your team.

- Clinical in-services
- Educational videos
- Tools, guides, checklists
- Downloadable information
- Ongoing clinical support



*\*Training kit is not intended for human use*



## EvoEndo® Training Kit\*

- EvoEndo® Head Model
- EvoEndo® Demo Controller
- Access to EvoEndo® online training platform

We achieve **big things** through **small scopes.**

Our ultra-slim technology unleashes the potential of sedation-free endoscopy to improve the lives of patients, families, and healthcare professionals.



## Contact Us for More Information

TEL 303.223.7445  
EMAIL [sales@evoendo.com](mailto:sales@evoendo.com)  
FAX 303.265.9381

[evoendo.com](http://evoendo.com)

  @evoendo