

Benefits

+ Safe

- The ONLY ultra-slim single-use gastroscope cleared by the FDA specifically for Upper GI endoscopy in Pediatric patients ages 5+.
- Enables Sedation-Free TNE, eliminating anesthesia risk.
- Single-use design minimizes cross-contamination.

+ Capable

- Available in 85 cm and 110 cm lengths for increased precision.
- 2.0 mm channel for standard pediatric accessories.
- 4-way steering with air, water, and suction functionality.

+ Flexible

- Portability makes system set-up a breeze.
- Low start-up costs make Sedation-Free TNE possible anywhere.
- Maximizes workflow efficiencies.



EvoEndo® Controller
Lightweight and portable
at 6" x 8.5" and 2 lbs.

Ready to learn more?

Scan the QR code to check out our website and schedule a hands-on demonstration.

On our website, you will find Clinical evidence, instructional videos, webinars, education opportunities, and more!



EvoEndo® Single-Use Endoscopy System v2.0



We achieve
big things
through
small scopes.

The EvoEndo® System is the only FDA-cleared endoscopy product designed **specifically for pediatric patients**, introducing a sedation-free, lower-risk, cost-effective alternative to traditional Upper GI endoscopy.

Our ultra-slim technology unleashes the potential of sedation-free endoscopy to improve the lives of patients, families, and healthcare professionals.



EMAIL info@evoendo.com

TEL 303.223.7445

FAX 303.265.9381

WEB evoendo.com

evoendo


evoendoinc

evoendo



Procedural Flexibility

Now with two scope lengths: a 110 cm and a softer and more flexible 85 cm option, taking the capabilities of the EvoEndo® Endoscopy System to the next level.




Enhance Image Feature

Unmatched definition and clarity in a single-use disposable scope.



Ergonomic Design

Familiar control layout fits a range of hand sizes for precision and confidence.



Optical System

Field of View	120°
Illumination method	LED

Insertion Portion

Bending section	210° Up, 90° Down, 100° Left, 100° Right
Maximum diameter of insertion portion	3.5 mm (0.14")
Working length	Model LE 85 = 85 cm Model LE 110 = 110 cm

Connections

Air Connector	Connects to 1/4"- 3/8" supply lines
Water Connector	Connects to bottle with 1.25" top
Suction Connector	Connects to 6 mm suction device

Channel

Average inner diameter	2.0 mm (0.078")
------------------------	-----------------

Model LE 85
85 cm

Model LE 110
110 cm



Patient-Centric

Our ultra-slim scope was created with small noses and narrow spaces in mind.

EvoEndo offers an Android-based virtual reality app* with single-use VR-goggles and a Patient Experience Kit. When used together, they provide patient distraction and can reduce anxiety.



*Requires Android OS mobile device with active wi-fi connection