



# PARASITOLOGY

# PARASITOLOGY

## Contents

- |   |                             |    |                           |
|---|-----------------------------|----|---------------------------|
| 3 | Collection and Transports   | 10 | Control Slides for Stains |
| 5 | Fecal Concentrating Systems | 11 | Parasite Suspensions      |
| 6 | Centrifuge Tubes            | 12 | Culture Media for Amoeba  |
| 7 | Stains and Reagents         | 13 | Microscope Supplies       |
| 9 | Staining Accessories        | 14 | Rapid Tests               |



Taenia solium egg



## At Hardy Diagnostics...

You're not just a number. You're not a figure on a graph in a quarterly report. At Hardy Diagnostics, you're a partner. From laboratories that utilize our tests to diagnose illness, to our employee owners who ensure that every lot meets your rigorous expectations, Hardy Diagnostics is about a group of people coming together to better the world, one test at a time.

Hardy Diagnostics is a 100% employee-owned company that is licensed by the FDA as a medical device manufacturer, and its quality management system is ISO 13485 certified. Hardy Diagnostics has been helping people live healthier lives since 1980.

Our microbiology products are used all over the world to diagnose and prevent disease. The dedication and skill of the people who work here have resulted in remarkable success.

Hardy Diagnostics donates 1% of its earnings to non-profit organizations and proudly supports a variety of charities and organizations in our communities.

In 2023, Hardy Diagnostics was awarded with the Premier GPO contract #PP-LA-640 and the U.S. Government GSA contract #47QSWA22D0077.



# Collection and Transports

## C&S Medium, Modified Cary Blair

Culture & Sensitivity medium provides a method for collecting and storing fecal specimens for the culture of intestinal enteric bacteria. Recommended for the recovery of *Shigella*, *Salmonella*, *Campylobacter*, and *Vibrio* species. A pH indicator is included to detect acidic conditions which are harmful for the recovery of enteric pathogens.

100/case 280505

100/case (with bags) 280505WB

## Formalin 10%

Neutral Buffered Formalin is an aqueous solution of formaldehyde buffered with sodium phosphate.

1 gallon, each 575A128

## PVA Fixative, Modified

Modified PolyVinyl Alcohol (PVA), an excellent alternative to mercuric chloride PVA for permanent staining.

32 oz, each VPM032

## ParaPak® Modified PVA/10% Formalin

Two vial system designed for the collection, transport, and preservation of stool specimens for intestinal parasites. Modified PVA Fixative for permanent staining that contains copper sulfate instead of mercuric chloride to eliminate the disposal problems associated with mercury. 10% Formalin (Neutral Buffered Formalin) is the recommended fixative for preserving fecal specimens for concentration and direct examination.

120 kits/case 300812

## ParaPak®, Ecofix®

This ecologically safer fixative contains no formaldehyde or mercury. It is suitable for direct examination, concentration, and permanent stain from the same specimen.

120/case 901312



### ParaPak®, Enteric Plus

This system is a modification of Cary Blair, designed for the collection, transport, preservation, and culture of bacterial enteric pathogens from fecal specimens.

120/case 900812

---

### Pinworm Paddle Collector

A clear polystyrene paddle with light adhesive on one side, allowing easy collection of pinworm and ova for microscopic examination. Features an integrated collector cap and crack-resistant polypropylene vial.

150/case V302PIN

---

### Total-Fix™/C&S /Clean Fixative Vial Set

The Total-Fix™/C&S/Clean fixative is a three vial system for the preservation of fecal specimens.

36 kits/case 89117

---

### Total-Fix™ Fixative

The Total-Fix™ fixative is a single vial system for the preservation of fecal specimens. Total-Fix™ preserved fecal specimens can be used for formalin sedimentation concentration, permanent stained smears, and most fecal immunoassays, such as DFA for *Cryptosporidium* spp. and *Giardia lamblia*.

100/case 280705

---

### Zinc PVA/Formalin Vial Set

An excellent alternative to mercuric chloride PVA for permanent staining. Formalin 10% fixative is recommended for preservation of fecal material used in concentration.

50 kits/case 89006

---

### Zinc PVA, Parasitology Fixative

Zinc Sulfate PVA-Mercury-free fixative. An excellent alternative to mercuric chloride PVA for permanent staining.

16 oz, each 2803A16

---



Helminths worms and their eggs

# Fecal Concentrating Systems

## ParaCon™ Filter Funnels

The ParaCon™ filter-funnel device provides a standard and consistent mesh size and allows filtration of fecal material without clogging.

96 funnels/per box PK96

1,000 funnels/per case PK1000



## Para-Pak® CON-Trate® System


A complete system for concentrating and recovering helminth eggs, larvae and protozoan cysts from feces.

Contains 50 filtering devices, 50 disposable centrifuge tubes with caps, Reagent A (MucoPenX) in a 15ml dropper bottle and Reagent B (Ethyl Acetate) in a 225ml container.

50 sets/case 960050



## Hookworm egg



# Centrifuge Tubes

## Sterile Centrifuge Tube

Blue screw cap, black printed graduations, white printed writing space, non-skirted, racked, sterile.

Material:

Tube: Polypropylene (PP)

Cap: High Density Polyethylene (HDPE)

Spin rate: 8,500 RPM

Temperature range: -20°C to +121°C

Centrifuge Tube, Sterile Racked,

15ml, 50 tubes/rack,

500/case 6286

50ml, 25 tubes/rack,

500/case 6289

15ml Centrifuge Tube



50ml Centrifuge Tube

## Non-Sterile Centrifuge Tube

Blue screw cap, black printed graduations, white printed writing space, non-skirted, bulk packed, non-sterile.

Material:

Tube: Polypropylene (PP)

Cap: High Density Polyethylene (HDPE)

Spin rate: 8,500 RPM

Temperature range: -20°C to +121°C

15ml, 500/case 6284

50ml, 500/case 6287



Cat. no. 6286

## Centrifuge Tube, Self-Standing

50ml centrifuge tube with separate red screw cap, skirted for self-standing, printed white graduation marks, non-sterile, bulk-packed.

Material:

Tube: Polypropylene (PP)

Cap: High Density Polyethylene (HDPE)

Spin rate: 5,000 RPM

Autoclavable (without screw cap)

50ml, 500/case 6255



# Stains & Reagents

## Acid Alcohol 90%

Acetic acid in alcohol solution. Decolorizer, a component of the trichrome procedure.

32 oz, each 3720A32

## Carbol Xylene

16 oz, each VXZ016

32 oz, each VXZ032

## Cryptosporidium Stain, Modified (Partial)

Partial acid-fast stain kit recommended for *Cryptosporidium* and *Isospora* species.

Partial AFB Stain Kit with Methylene Blue counterstain  
4 x 8 oz bottles/kit 482K

Partial AFB Stain Kit with Malachite Green counterstain  
4 x 8 oz bottles/kit 482KMG

## Ethanol 70%

For trichrome and other uses, reagent grade.

32 oz, each VXB032

One gallon, each 3716A128

## Envirene, Xylene Substitute

An odorless (no citrus smell) non-toxic substitute for xylene or procedures that call for xylene, such as histology and parasitology.

16 oz, each CE016

1 gallon, each CE128

## Ethyl Acetate

For preparing parasite concentrations.  
16 oz, each ET105

## Ethanol 95%

For trichrome stains and other uses.  
32 oz, each VXC032

1 gallon, each 3719A128

## Iodine Alcohol

For removing mercury deposits in Schaudinn or PVA fixed fecal smears.

16 oz, each VXI016

1 gallon, each VXI128

Kinyoun Carbol Fuchsin with Methylene Blue counterstain, Cat. no. 482K





Wheatley's Trichrome Stain  
16 oz, Cat. no. T016

### Schaudinn's Fixative, Modified

Mercury-free (copper sulfate formula) and used as an on-slide parasite fixative for fresh fecal smears.

Modified Schaudinn's Fixative is used in Wheatly trichrome stain procedures and prevents distortion of parasitic elements.

16 oz, each VXM016

32 oz, each VXM032

### Trichrome Starter Kit

Trichrome parasitology starter kit (Wheatley's Procedure) has everything you need to get going. This single, convenient package contains all the materials to prepare stained fecal smears.

The kit contains:

- 8 ounces each of Schaudinn's Fixative and Trichrome Stain
- 16 ounces each of Carbol Xylene, Xylene, Iodine Alcohol, 70% Ethanol, 95% Ethanol, 90% Acid Ethanol
- 10 reusable Coplin Jars

Additionally, Hardy sells a refill kit (Cat. no. VXR010) with all of the above listed items except Coplin Jars.

Kit VXF019

### Iodine, Dobell and O'Connor

A weak iodine solution recommended for use in staining protozoan cysts. This preparation is a 1:5 dilution of Lugol's iodine.

Amber glass bottle, 15ml Z66

### Mayer's Albumin

Used to adhere fecal specimens to glass slides in the Modified Iron Hematoxylin/Kinyoun or trichrome staining process.

15ml, each Z69

### Modified Trichrome Blue Stain

For staining microsporidia, isospora, and cyclospora. Contains aniline blue.

16 oz, each MT016

# Staining Accessories

## Coplin Jar, Plastic

For manual staining procedures, jars are white plastic with a screw-on top to minimize evaporation. These jars are for use in all staining procedures, with slots on top to eliminate slide contact. Each jar holds 5 slides in upright position or 10 slides back-to-back.


Each VCJ001



## Slide Staining Rack

This stainless steel staining rack is designed for use with laboratory sinks of various sizes. Triangular slide supporting rods have a tapered top and contact the slides only in a narrow zone. Rack can be easily and quickly adjusted to fit sinks up to 53.4cm (21") inside dimensions.

24" long, each 15153821



# Control Slides for Stains

## Acid-Fast Control Slide

This quality control slide provides two air-dried and methanol-fixed droplets within two etched circles. The circle nearest the label contains a droplet of an acid-fast-positive *Mycobacterium*. The circle furthest from the label contains a droplet of an acid-fast-positive *Cryptosporidium* in a fecal smear sample and contains acid-fast-negative intestinal bacteria.

10 slides/box SL4010

## Blood Parasite Control Slide

This quality control slide provides an air-dried, methanol-fixed, blood parasite smear containing *Plasmodium*, *Babesia*, or *Trypanosoma*.

10 slides/box SL9010

## Cryptosporidium Control Slide

This quality control slide provides a single air-dried and methanol-fixed fecal smear containing acid-fast-positive *Cryptosporidium* and acid-fast-negative intestinal bacteria.

10 slides/box SL4510

## Intestinal Protozoa, SAF, Control Slide

Each slide has a single air-dried, SAF-preserved fecal smear with a representative intestinal protozoa, usually *Giardia lamblia*.

10 slides/box SL1510

## Intestinal Protozoa, Zinc PVA Control Slide

Each slide has a single air-dried, Zinc PVA preserved fecal smear with a representative intestinal protozoa, usually *Giardia lamblia*.

10 slides/box SL1010

## *Pneumocystis carinii* Control Slide

This quality control slide provides a single air-dried and methanol-fixed non-human tissue preparation containing *Pneumocystis*.

10 slides/box SL7010

## Q-Slide™ AFB Control Slide

A convenient ready-to-stain control slide that will test both your stain solutions and your staining techniques. Positive and negative stain control organisms are affixed to the circles in each Q-Slide™.

5/pk Z301

Acid-Fast Control Slide,  
Cat. no. SL4010



# Parasite Suspensions for QC

Parasite suspensions are packaged in a vial containing 1ml of a particular parasite species. A Certificate of Analysis is available for each lot of parasite suspension that includes detailed information about the parasite. Choose from these common microorganisms below:

<b>Parasitology</b>	FP01	<i>Giardia lamblia</i> Parasite Suspension
	FP02	<i>Cryptosporidium</i> species Parasite Suspension
	FP03	<i>Diphyllobothrium latum</i> Parasite Suspension
	FP04	<i>Taenia</i> species Parasite Suspension
	FP05	<i>Ascaris lumbricoides</i> Parasite Suspension
	FP06	<i>Strongyloides stercoralis</i> Parasite Suspension
	FP08	<i>Trichuris trichiura</i> Parasite Suspension
	FP09	<i>Iodamoeba buetschlii</i> Parasite Suspension
	FP10	<i>Necator americanus</i> Parasite Suspension
	FP11	<i>Hymenolepis nana</i> Parasite Suspension
	FP12	<i>Endolimax nana</i> Parasite Suspension



## Features

- Ready-to-use
- Hassle-free packaging
- Room temperature storage
- Certificate of Analysis
- CE Marked medical IVD

## Applications

- QC of diagnostics kits and methods
- Microscopic Examination
- Traditional staining procedures
- Proficiency programs
- Education

# Culture Media for Amoeba

## Free-Living Amoeba Media

For the detection of *Acanthamoeba* and *Naegleria* species from clinical and environmental samples. Microscopic examination beginning on the second day of incubation reveals the presence of typical “feeding tracks,” which indicates presence of amoebae. Requires *E. coli* stock culture (not provided) for use in preparing the bacterial lawn on the Non-Nutrient Agar plates.

## Acanthamoeba Broth (PYG)

For the maintenance of *Acanthamoeba* stock cultures.

16x100mm tube, 5ml,  
20/pk

K225

## Page's Saline

For the preparation of suspension used to inoculate Non-Nutrient agar plates. 13x100mm polycarbonate tube, 2ml, 20/pk


R225

## Non-Nutrient Agar

Non-Nutrient Agar plate for the detection of free-living amoebae including *Acanthamoeba* and *Naegleria* spp.

15x100mm plate, 10/pk

G225



*Acanthamoeba castellanii* showing “feeding tracks” on Non-Nutrient Agar with *E. coli*.

Non-Nutrient Agar,  
Cat. no. G225

# Microscope Supplies

## Bibulous Blotting Paper

High quality, highly absorbent papers manufactured for minimum linting and shedding. Recommended for drying slides and small parts or components. 50 sheets per book.

12 books/pack 28511007

## Cover Glass

High-quality, borosilicate cover glass is non-corrosive and non-fogging and is devoid of scratches, bubbles, and striations. The flatness provides a uniform distortion-free surface. A desiccant is added to each box to ensure that moisture doesn't degrade the cover glass. Convenient plastic storage container allows for easy stacking on the shelf.

18x18mm, 1 oz box 1818HD

22x22mm, 1 oz box 2222HD

22x40mm, 1 oz box 2240HD

24x40mm, 1 oz box 2440HD

24x50mm, 1 oz box 2450HD

## ProSlide™ Microscope Slides

Single frosted one end on one side. White glass, beveled edges. Available in the super clean and polished premium grade or economy grade. The economy slide is not recommended for hematology use.

3x1 inches x 1mm, 72/box:

frosted, economy PF72E

frosted, premium PF72P

plain, economy PP72E

plain, premium PP72P

## Lens Cleaner

An effective microscope lens cleaning solution in a dropper bottle. Safe to use on the most delicate optics. Use with lens tissue, Cat. no. 52846001.

30ml, each Z97

## Bibulous Blotting Paper



## Lens Paper, Non-linting

For cleaning microscope lenses, safety glasses, goggles, and more with no scratching. Dust-free lens paper tissues are ideal for use with liquid cleaners.

4x6 inches, 50 sheets/pk,

12/pks 52846001

## Immersion Oil

Used to assist in microscopic examination at high power by increasing the resolving power.

### Low Viscosity

Type A, 20ml, each Z95

Type A, 4 oz, each Z85

### High Viscosity

Type B, 20ml, each Z96

Type B, 4 oz, each Z86

## Immersion Oil





## Rapid Tests

### PARA-TECT™ Cryptosporidium/ Giardia Direct Fluorescent Assay

This direct fluorescent assay (DFA) is an *in vitro* immunoassay for the qualitative determination of *Cryptosporidium* and *Giardia* oocysts and cysts in feces. Can be used with stools preserved in 10% formalin, SAF, or TOTAL-FIX™ Universal Fixative.  
75 tests/pack MCCCG75

### Merifluor®, Crypto/Giardia

Merifluor® is an *in vitro* direct immunofluorescent detection procedure for the simultaneous detection of *Cryptosporidium* oocysts and *Giardia* cysts in fecal material. Requires just 30 minutes of incubation. Bright apple-green fluorescence of the parasite against red-orange background makes rare microorganisms easy to see.  
50 tests/pack 250050





# HARDY DIAGNOSTICS

*A Culture of Service™*



40 Years of  
Bringing You  
A Culture of Service



100%  
Employee Owned



Hardy Diagnostics  
donates 1% of each  
sale to charity.



FM 572526

Hardy Diagnostics has a Quality Management System that is certified to ISO 13485 and is a FDA licensed medical device manufacturer.

## Headquarters

1430 West McCoy Lane  
Santa Maria, CA 93455  
800.266.2222  
Sales@HardyDiagnostics.com  
HardyDiagnostics.com

## Distribution Centers

Santa Maria, California  
Olympia, Washington  
Salt Lake City, Utah  
Phoenix, Arizona  
Dallas, Texas  
Springboro, Ohio  
Lake City, Florida  
Albany, New York  
Raleigh, North Carolina



# PREMIER®

Contract # PP-LA-640 | **Contracted Supplier**



# Contract Holder

Contract # 47QSWA22D0077