



Hygiena's Digital Dry Block Incubators are designed specifically for Hygiena's wide variety of test devices. Capable of up to 105 °C, the incubators can be set to any temperature required for Hygiena products. The digitally controlled, interchangeable blocks of the Lab Format Incubator are independently heated, allowing you to run two temperatures on one incubator at the same time – ideal for running a variety of tests. Auto-stop and sound alert features help you achieve perfect incubation time, every time. Available in two sizes and up to 70 wells, Hygiena's Digital Dry Block incubator is configurable for your unique needs.


Benefits:

- Designed specifically for Hygiena test devices
- Numbered wells help keep testing organized
- Flexible configuration can accommodate swab test devices or large volume tubes
- Digitally controlled independent blocks can be held at different temperatures
- Small footprint takes minimal lab bench space

Features:

- Wide temperature range:
 - RT+5 105 °C for Lab Format / RT+5 80 °C for Small Format Incubator
- Temperature stability $\leq \pm 0.5$ °C
- Heats to max temp in less than 15 minutes
- Timer range: 1 min 99 hr 59 min
- Footprint / Weight:
 - Small Format (WxDxH): $11 \text{ cm} \times 15 \text{ cm} \times 8 \text{ cm}$, $0.5 \text{ kg} / 4.3 \text{ in} \times 5.9 \text{ in} \times 3.1 \text{ in}$, 1.1 lbs
 - Lab Format (WxDxH): $36.5 \text{ cm} \times 21 \text{ cm} \times 15 \text{ cm}$, $4.5 \text{ kg} / 14.4 \text{ in} \times 8.3 \text{ in} \times 5.9 \text{ in}$, 9.9 lbs















Available Configurations for Lab Format Incubator:




Lab Format Incubator

Catalog No.	Description	Qty
INCUBATOR2	Lab Format Digital Dry Block Incubator	1
IB001	Hygiena Swab Tube Block – 35 wells	1
IB002	Hygiena Vial Tube Block – 15 wells	1 1







(1) IB001 & (1) IB002

(2) IB002

Available Configuration for Small Format Incubator:



Small Format Incubator

Catalog No .	Description	Qty
INCUBATOR	Small Format Digital Dry Block Incubator	1
IB003	Hygiena Swab Tube Block – 12 wells	1



(1) IB003



Headquarters

1430 West McCoy Lane Santa Maria, CA 93455 800.266.2222 :phone 805.346.2760 :fax Sales@Hardy Diagnostics.comwww.HardyDiagnostics.com

Distribution Centers

REV A 052014

Santa Maria, California Olympia, Washington Salt Lake City, Utah Phoenix, Arizona Dallas, Texas

Springboro, Ohio Lake City, Florida Raleigh, N. Carolina Albany, New York