
Hardy CHROM[™] Group A Strep Agar






Hardy**CHROM**[™] Group A Strep Agar

Color read-out on a white background in just 24 hours.

We are making identification easy with the release of the Hardy Diagnostics HardyCHROM[™] Group A Strep Agar. This easy-to-read chromogenic medium is recommended for the selective cultivation and differential isolation of Group A *Streptococcus* (*S. pyogenes*). Colonies are identified based on colony color (red, red-brown, or red-orange colonies) among the other non-GAS bacteria in the complex throat flora (blue, clear or white colonies) after 24 hours of incubation.

15x100mm plate, 10/pk	G337
Reduced stacking ring,	
15x100mm plate, 10/pk	GA337




GROUP A STREP infections can occur any time during the year. However, infections are more common in the United States in certain seasons:

- Strep throat and scarlet fever are more common in the winter and spring.
- Impetigo is more common in the summer.

Anyone can get a **GROUP A STREP** infection, but some infections are more common in certain age groups:

- Strep throat and scarlet fever are most common in children between the ages of 5 and 15 years.
- Impetigo is most common in children between the ages of 2 and 5 years. (source CDC)



Streptococcus pyogenes (ATCC® 19615) red-brown colonies growing on HardyCHROM™ Group A Strep Agar. Incubated in CO₂ for 24 hours at 35°C.



Streptococcus pyogenes (clinical strain) red-orange colonies growing on HardyCHROM™ Group A Strep Agar. Incubated in CO₂ for 24 hours at 35°C.

Did you know that Hardy Diagnostics offers you free blok-boxes for all your light sensitive media? Just ask us if you would like a free blok-box added to your order. One box holds a sleeve of ten Petri plates.









HUSA





Hardy Diag<mark>nostics</mark> donates 1% of each sale to charity.



Hardy Diagnostics has a Quality Management System that is certified to ISO 13485 and is a FDA licensed medical device manufacturer.

Headquarters

1430 West McCoy Lane Santa Maria, CA 93455 800.266.2222 Sales@HardyDiagnostics.com HardyDiagnostics.com

Distribution Centers

Santa Maria, California Olympia, Washington Salt Lake City, Utah Phoenix, Arizona Dallas, Texas Springboro, Ohio Lake City, Florida Albany, New York Raleigh, North Carolina

