Infection control and patient care.

They're one and the same.





Witnessefficient*Clostridioidesdifficile* detection in as little as 24 hours!



Hardy Diagnostics introduces C diff Banana Broth[™], which is recommended as a selective detection broth for the recovery of *Clostridioides difficile* vegetative cells or spores from environmental samples. 20/pk K226

- Specificity found to be 100%. Ask for reference
- No special anaerobic incubation equipment required
- Selective and differential medium for recovering *Clostridioides difficile* spores and vegetative cells
- No need to pre-reduce the media
- Easy color change read-out





Negative

Positive

Verifytheefficacyofyourcleaning protocol




Saline

For pre-moistening of the swab. 0.85%, 2ml, 13x100mm tube, shatter resistant, polycarbonate tube. 20/pk R45

UNE 0.85% R45 13046 2018-02 RDY

Environmentaldetection of Clostridioides difficile



KWIK-STIK[™] (Positive Control)

A ready-to-use swab that contains a single lyophilized pellet of a single microorganism strain. *Clostridioides difficile*, ATCC^{*} 9689^{°°}. To use, simply crush the rehydrating ampule, mix, remove swab, and inoculate, by Microbiologics^{*}. 6/pk 0329K



Boekel Incubator

- Dimensions: 11.62" x 15.5" x 15.5"
- Temp. range: Ambient +10°C to 90°C
- LCD temperature indicator
- Supplied with a single removable zinc plated steel shelf

Each

helf 138000













Hardy Diagnostics donates 1% of each sale to charity.



FM 572526 agnostics has a G

Hardy Diagnostics has a Quality Management System that is certified to ISO 13485 and is a FDA licensed medical device manufacturer.

Headquarters

1430 West McCoy Lane Santa Maria, CA 93455 800.266.2222 Sales@HardyDiagnostics.com HardyDiagnostics.com

Distribution Centers

Santa Maria, California Olympia, Washington Salt Lake City, Utah Phoenix, Arizona Dallas, Texas Springboro, Ohio Lake City, Florida Albany, New York Raleigh, North Carolina

