DAILY SHIFT HEADS | For use with TRIO.BAS™

Suitable for:

- Clean Room
- Isolator/RABS
- HVAC
- Field Testing
- Indoor Air Quality
- · Air Filter Monitoring

For use with either 100mm or 60mm plate format TRIO.BAS™ instruments:

- MINI
- MONO
- DUO
- **ISOLATOR**

Daily Shift Head Features:

- Design is identical to stainless steel counterparts
- · Sterile: gamma irradiated to assurance level SAL 109
- Certificateofsterility included
- · Individual Daily Shift Heads are double bagged
- · Labeled with lot number and expiration date
- 5 Year expiration from date of manufacture
- Constructed of antistatic resin plastic
- 100% Recyclable
- Packed in a box of 27





Cat. no.

DAILY SHIFT HEAD, Contact plate, 20/bx

BAS340

DAILY SHIFT HEAD, Petri plate, 20/bx

BAS341





Why Adopt Daily Shift Heads?:

- Daily Shift Heads can be used on the same day, on the same shift, inside the same cleanroom
- The use of a sterile aspirating head with an active air sampler is requested for each microbial air sampling cycle, or group of sampling cycles, in the same controlled environment
- Stainless Steel Aspirating heads should be sterilized by autoclaving with subsequent filing of a sterilization document. Daily Shift Heads are ready-to-use.
- Limit possible operator contamination
- Save operator time and sample with better efficiency, avoiding cleaning, autoclaving and validation
- Have on hand to easily present a certificate of sterility assurance to a regulatory authority
- Have on hand in the case of autoclave failure so sampling is not disrupted
- Have on hand when sampling demands are especially high







Contact your Sales Representative or visit http://hardydiagnostics.com/triobas/
800.266.2222 Sales@HardyDiagnostics.com