# Products to Control COVID-19 and INFLUENZA





#### Contents

- 1 Transport Medium
- 2 Rapid & LFA Tests
- 3 Influenza Test Kit
- 4 Swabs
- 4 Equipment
- 5 Disinfectants & Sanitizers
- 5 Personal Protective Equipment

DIAGNOSTICS A Culture of Service<sup>™</sup>

Be sure to visit our website at HardyDiagnostics.com/coronavirus-covid-19-update/ for all the lastest news and products.



# **Transport Medium**

# TransPRO VTM

### Kit For all your viral transport needs. From Hardy

#### Diagnostics, the partner you know and trust!

TransPRO<sup>™</sup> VTM is recommended for the preservation of viral agents including Influenza A, Influenza B, Adenovirus, and Echovirus.

The transport medium comes in a plastic, screw cap tube with self-standing skirt and inner conical base.

- Manufactured in the USA
- Room temperature\* stable
- Each kit contains one (1) tube of TransPRO<sup>™</sup> VTM (3ml fill) and one (1) individually wrapped (sterile) nasopharyngeal swab. All packaged in an easy to use specimen collection peel pouch kit
- Viral transport media contains antibiotics and antifungal agents to inhibit bacterial and fungal growth within patient samples

\* See IFU for further details 50/pk





## Viral Transport Medium (VTM)

**Intended use**: For the collection and transport of clinical specimens for the preservation of viral agents including, Influenza A, Influenza B, Adenovirus, and Echovirus from the collection site to the testing laboratory. Hardy Diagnostics' VTM is a culture-based media that is intended to be used in standard laboratory procedures for virus culture and diagnostic assays that utilize stable recoverable infectious viral particles.

**Storage:** Internal studies support that the media can be stored at refrigerated or room temperature storage (2-30°C) for up to 365 days without decreased performance. Data permitting, the length of RT storage shall continue to be extended as the realtime findings are collected, and timely updated in the product IFU.

17x100mm polypropylene tube, 3ml fill, 20/pk

TRV50

R99



# **Rapid & LFA Tests**

### COVID-19 IgG/IgM Rapid Antibody Test bv RightSign<sup>™</sup>

The RightSign COVID-19 IgG/IgM Rapid Test cassette is a lateral flow immunochromatographic assav for the detection of SARS CoV-2 antibodies in venous whole blood, serum, plasma, and fingerstick whole blood. The COVID-19 IgG/IgM Rapid Test is intended for use as an aid in identifying individuals with an adaptive immune response (due to infection or vaccination) to SARS-CoV-2, indicating recent or prior infection.

#### **CLIA WAIVED**



#### Features and Benefits

- Identifies IgG and IgM antibodies to SARS-CoV-2 with a combined sensitivity and specificity of 100%.
- Read results in only 10 minutes.
- No equipment needed.

**EUA Authorized** 

Sensitivity and Specificity

Visually read in 15 minutes

Specificity 100%

Type B Antigen.

25 tests/bx

- Can use fingerstick whole blood for Point of Care (POC) testing (CLIA Waived).
- IgM indicates a current COVID-19 infection and IgG indicates a past infection.

20/pk

**RTCV1920** 



#### Status<sup>™</sup> Flu A&B Test Kit, CLIA Waived

The LifeSign Status Flu A & B Test Kit is used by healthcare professionals to provide an easy and convenient point-of-care testing of Influenza A or B virus in a symptomatic patient. Meets the new FDA reclassification criteria and CLIA-waived, specimen needed for testing can be collected directly through anterior nasal or nasopharyngeal swabbing. Accurate and fast results will be shown on the test device in as fast as 10-15 minutes, and promotes efficiency and precision allowing proper differentiation of Flu A and B infection.

Safe & Reliable Kit: Each in-vitro rapid testincludes a nasal swab, developer solutioncapsule, and a test cassette, providing simpleand easy testing. The flock nasal swab enhancesthe specimen collection, while the pre-measureddeveloper solution capsule improves testaccuracy and ease of use.25 tests/bx36025

#### Status<sup>™</sup> Flu A&B Control Swabs

The LifeSign Status Flu A & B Positive and Negative Control Swabs are external controls that guarantee medical professionals optimal performance of the Status Flu A & B Test Kit.

Ensures Accuracy and Reliability. Designed to check the performance of the Test Kit, the Status Flu A & B Positive and Negative Control Swabs contain recombinant protein and inactivated microorganisms. The Positive Control Swab has recombinant nucleoproteins of influenza A and B viruses that are dried onto the swab while the Negative Control Swab consists of inactivated Group B Streptococcus dried onto the swab.

25 swabs/pk

36022CS







# Influenza Test Kit

### QuickVue<sup>™</sup> Influenza A+B Kit

Tests for and differentiates between Influenza A or B, 30 seconds of hands-on time, 10 minutes total time, and easy readout. CLIA waived.

25 tests/kit





#### Flocked Nasopharangeal (NP) Swab, Sterile, 80mm Breakpoint

Flexible, mini-tip size flocked swab with 80mm breakpoint in peel pouch, individually packaged, sterile, Nasopharyngeal (NP).



#### Mid-Turbinate Adult flocked swab w/stopper, 80mm Breakpoint

Flocked tip swab, 80mm breakpoint, in peel pouch - individually packaged, sterile, disposable.



#### Mid-Turbinate Pediatric flocked swab w/stopper, 80mm Breakpoint

Regular polyester tip swab, 80mm breakpoint, in peel pouch - individually packaged, sterile, disposable.





## Equipment

#### UV Chamber, N95 Mask Decontamination by UV light

Never run out of N95 masks again. Quickly sanitize them yourself by UV light.

- Sanitizes four N95 masks at once. Takes only two minutes.
- Two modes of operation Energy or Time
- Open door shut-off safety mechanism
- Energy efficient electronic ballasts
- Small footprint
- Two Year Manufacturer's Warranty
- UL/CSA/CE Approved
- Equipped with 5 x 254 nm Shortwave Bulbs



234100HD

Each

# **Disinfectants & Sanitizers**

#### **Alpet D2 Surface** Sanitizer/Disinfectant

- 58% isopropyl + quaternary ammonia
- Clean, sanitize and disinfect surfaces with one product
- 60 second kill on food contact surfaces, including gloves
- 10 second kill on non-food contact surfaces
- Highly evaporative and ideal for water sensitive equipment and low moisture environments



SS1000312

#### 1L bottle. 12/pk

#### **Bio-Pure**<sup>™</sup> Alcohol-Wipes

- 70% IPA /30% H<sub>2</sub>O mixture evaporates completely leaving no residue
- Wipes won't scratch surfaces of optics
- Lint free
- Strong, durable fabric
- 100% biodegradable



#### 40 wipes, Each

**BPWP1000** 

#### **CaviWipes**<sup>™</sup> Disinfectant **Towelettes**

- Presaturated with CaviCide<sup>®</sup>
- CaviWipes<sup>™</sup> are the only cloths proven to kill both TB (in 5 minutes) and HBV (in 3 minutes)
- Quick, easy-to-use
- Timesaving convenience

160 wipes. Fach

#### AloeSafe<sup>™</sup> Antiseptic Gel Hand Sanitizer

AloeSafe<sup>™</sup> hand sanitize with ethyl alcohol gel with Aloe Vera, Lavender Oil, and Vitamin E, without the harsh chemicals.

- Formulated with Aloe and Vitamin E to keep from drying
- Lavender Scent
- 62% Ethanol
- No Glycerin
- No Parabens
- Kill 99.9% of Bacteria

16 oz. bottle. Each



131100

### AloeSate Antiseptic Hand Sanitize HARDY

AS16HD

Personal Protective Equipment

#### Adenna Disposable Face Mask

4-Fold Fluid Resistant Earloop Face Masks offer a dependable solution in applications whereby maximum coverage is preferred. These masks

are 100% latex and fiberglass free, as well as non-irritant and odorless. The unique consealed PVC nose piece conforms to your face better for a more comfortable wear. 10x50 case/500 masks



### Face Shield, Full Face, Disposable

Face shields are personal protective equipment (PPE) that are used for protection of the facial areas (eves, nose, mouth) from sprays, aerosols. and splashes.

- Made from PETG.
- Provides clarity of view
- Anti-fogging visor
- Once size fits most
- Comfortable Latex Free forehead cushion

10/pk

MD30331















Hardy Diagnostics donates 1% of net profits to charity.



Hardy Diagnostics has a Quality Management System that is certified to ISO 13485 and is a FDA licensed medical device manufacturer.

#### Headquarters

1430 West McCoy Lane Santa Maria, CA 93455 800.266.2222 Sales@HardyDiagnostics.com HardyDiagnostics.com

#### **Distribution Centers**

Santa Maria, California Olympia, Washington Salt Lake City, Utah Phoenix, Arizona Dallas, Texas Des Moines, Iowa Springboro, Ohio Lake City, Florida Albany, New York Raleigh, North Carolina