


AlerTox[®] Sticks

Rapid tests for allergen detection in food, drinks, and surfaces for internal Quality Control.



Egg



Crustacean



Hazelnut



β-Lactoglobulin



Fish



Almond



Casein



Mustard Seed



Walnut



Total Milk



Soy



Peanut

AlerTox[®] Sticks advantages:

- » Quick and Reliable
- » Results in 10 minutes
- » No Cross-Reactivity
- » High Sensitivity

AlerTox[®] Sticks are a practical tool for the control of allergen content in raw materials, final products, and working surfaces.

Industrial application

- » Semi-quantification of allergens in food
- » Food quality control
- » Trace allergen contamination in food, liquid and surfaces
- » Implementation and fulfillment of HACCP, ISO 22000, BRC, and IFS compliance



Product Name	Cat. no.	Quantity
AlerTox [®] Sticks Egg	KIT3025	10 tests
AlerTox [®] Sticks Crustacean	KIT3036	10 tests
AlerTox [®] Sticks Hazelnut	KIT3035	10 tests
AlerTox [®] Sticks Fish	KIT3038	10 tests
AlerTox [®] Sticks β -Lactoglobulin	KIT3018	10 tests
AlerTox [®] Sticks Almond	KIT3033	10 tests
AlerTox [®] Sticks Mustard Seed	KIT3040	10 tests
AlerTox [®] Sticks Casein	KIT3021	10 tests
AlerTox [®] Sticks Total Milk	KIT3015	10 tests

Learn more at HardyDiagnostics.com/Food

NOTE: As with any qualitative or semiquantitative test, results should be periodically validated with a quantitative test performed in an outside laboratory.