Status Strep A

Rapid immunochromatographic assay for the detection of group A streptococcal antigen directly from throat swab specimens

Status Strep A-Flip Cassette



■ CLIA WAIVED

- Patented single use reagent capsules
- Positive results in 5 minutes
- Sensitivity of 96% and specificity of 99%
- No counting reagent drops -Less steps for greater accuracy
- External negative and positive controls included

CPT CODE 87880OW



Status Strep A Strip



CLIA WAIVED

- Simple procedure
- Easy extraction for fast reaction times
- Results in 5-10 minutes
- Sensitivity of 96.2% and specificity of 98.7%
- External negative and positive controls included

CPT CODE 87880QW Two convenient resealable canisters holding 25 test strips each

Status Strep A Plus

Strip





■ Non-waived

- Simple procedure
- Easy to read results
- Results in 5-10 minutes
- Sensitivity of 96.2% and specificity of 98.7%
- External negative and positive controls included

CPT CODE 87430

ORDERING INFORMATION			
PRODUCT	CLIA COMPLEXITY	# OF TESTS	ITEM #
Status Strep A Plus (Cannister-Strip	s) Waived	50	34250
Status Strep A (Strip)	Waived	30	34130
Status Strep A Flip (Cassette)	Waived	25	34125
Status AccuStrep A (Cassette)	Moderate	25	34025

Proudly made in the USA. For more information or to place an order, contact



Scan to find your local Hardy Sales Rep



Official Distributor of

